INSTRUCTIONS FOR USE

INION BIORESTORE™

DESCRIPTION

The **INION BIORESTORE™** is an osteoconductive and osteostimulative bioactive bone graft device that can be used either separately or in conjunction with autogenous or allograft bone*. In vivo tests have demonstrated more bone formation at each post-implantation time point and more total bone formation compared to other osteoconductive devices such as hydroxyapatite. In vitro cell culture tests with human adipose mesenchymal stem cells have demonstrated an osteostimulative effect, defined as the active stimulation of osteoblast proliferation and differentiation as evidenced by alkaline phosphatase activity. This stimulation has been attributed as being the result of the interaction between osteoblasts and the ionic dissolution products released from **INION BIORESTORE™** during its absorption. Clinical data in humans on rate and extent of bone formation observed in cell culture and animal models has not been established.

*Mixing with autogenous or allograft bone is not cleared in the USA for use in the extremities, spine and pelvis!

Bioactive materials are those materials that elicit a specific biological response at the interface of the material that results in the formation of a bond between the tissues and the material. Osteostimulation is a "...property of some bioactive materials to enhance, actively stimulate both the proliferation and differentiation of progenitor cells (e.g., mesenchymal stem cells)..."^{1,2}

¹ Society for Biomaterials presentation, 24th Annual Meeting, 1998:511-518

² J Appl Biomat 1992, 3:123-129

The **INION BIORESTORE™** consists of different size/shape porous blocks and morsels made of degradable bioactive glass. Bioactive glasses have a long history of safe medical use and undergo a time-dependent, kinetic modification of the surface that occurs when implanted in living tissue. Specifically, the surface reaction results in the formation of a calcium phosphate layer that is substantially equivalent in composition and structure to the hydroxyapatite found in bone mineral. This apatite layer provides scaffolding onto which the patient's new bone will grow allowing complete repair of the defect. Based on preclinical testing, most of the material degrades in vivo in six months. The material is radiopaque.

The **INION BIORESTORE™** implants are sterile and non-collagenous.

INDICATIONS

Indications in the EU

The **INION BIORESTORE™** implants are bone graft substitutes indicated for bony voids or gaps that are not intrinsic to the stability of the bony structure; packing into bony voids or gaps of the skeletal system to fill and/or augment dental intraosseous, oral and maxillofacial defects, and defects in the extremities, spine and pelvis. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone, including: e.g., periodontal/infrabony defects; alveolar ridge augmentation (sinusotomy, osteotomy, cystectomy); dental extraction sites (ridge maintenance, implant preparation/placement); sinus lifts; cystic defects; maxillofacial augmentation; acute long bone fractures; bony defects existing in either acetabulum, femur, or tibia; and defects associated with spinal stabilization and fusion.

Indications in the USA

BRS-1600 - 1620, and BRS-5400 - 5450:

These **INION BIORESTORE™** implants are bone graft substitutes indicated only for bony voids or gaps that are not intrinsic to the stability of the bony structure. **INION BIORESTORE™** is indicated to be gently packed into bony voids or gaps of the skeletal system (i.e., the extremities, spine and pelvis). These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone. The product provides a bone void filler that resorbs and is replaced with bone during the healing process.

BRS-1400 - 1410, and BRS-2100 - 2120:

These **INION BIORESTORE™** implants are bone graft substitutes indicated for bony voids or gaps that are not intrinsic to the stability of the bony structure; packing into bony voids or gaps to fill and/or augment dental intraosseous, oral and craniomaxillofacial defects. These defects may be surgically created osseous defects or osseous defects created from traumatic injury to the bone, including: e.g., periodontal/intrabony defects; alveolar ridge augmentation (sinusotomy, osteotomy, cystectomy); dental extraction sites (ridge maintenance, implant preparation/placement); sinus lifts; cystic defects; and craniofacial augmentation.

CONTRAINDICATIONS

The **INION BIORESTORE™** implants are not intended for use in and are contraindicated for:

- Active or potential infection.
- Patient conditions including:
 1. Use of medication known to affect the skeleton (e.g., chronic glucocorticoid usage >10mg/day for the previous 3 months). Estrogen replacement therapy is allowed.
 2. Need for chronic anticoagulant therapy (e.g., heparin). Prophylactic use of Coumadin or aspirin postoperatively is allowed.
 3. A systemic metabolic disorder known to adversely affect bone healing and mineralization (e.g., poorly controlled insulin-dependent diabetes, renal osteodystrophy, Paget's disease), other than primary osteoporosis.
 4. Any existing condition or disease that will interfere with good soft tissue and bone healing.
 5. Defects requiring graft material to exceed the volume of 30 cm³.

6. Limited blood supply.
7. When patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse).

INFORMATION FOR USE

Surgical considerations and reminders

- Prophylactic preoperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.
- Maintain sterile field throughout the procedure.
- Proper exposure using standard surgical procedure.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Make every attempt to obtain complete tension-free primary soft tissue closure over the material.

Preparation of the INION BIORESTORE™

- Choose appropriate INION BIORESTORE™ block(s) or morsels for the indication. Preoperative radiographic evaluation of the defect site is essential to accurately assess the extent of the defect and to aid in the selection and placement of the graft material.
- The INION BIORESTORE™ blocks may be cut or shaped to the desired three dimensional shape with surgical instruments (e.g., surgical knives, scalpels or scissors). However, before shaping the INION BIORESTORE™, the implant should be moistened with either sterile saline or water. Furthermore, the shaped implant should always be rinsed with sterile saline or water after shaping and before any contact with the patient's tissues to remove excessive loose particles from the surfaces of the device. Take care to ensure that any shaped device surfaces are smooth and free from excessive loose particles.
- The INION BIORESTORE™ should be wet with patient's blood or marrow, or sterile saline or water before application. The osseous regeneration may occur more rapidly by using the patient's osteogenic blood or marrow obtained by scoring the surface or drilling into the cortical bone at the surgical site. Place the desired fluid in a sterile basin or other suitable container, or inject blood directly into the syringe containing BioRestore morsels. Place the material in contact with the fluid to soak up the fluid. The porous structure of the material will allow the fluid to permeate into the implant by capillary action.

Implantation of the INION BIORESTORE™

- SITE PREPARATION:
 - Once the defect has been exposed, thorough debridement of the osseous defect of all necrotic and soft tissue should be carried out. Eliminate all granulation tissue or cysts at the defect site. The walls of the defect should be decorticated.
 - Irrigate the defect with sterile saline or water and evacuate the excess.
 - Some bleeding should be observed originating from the host bone to indicate viability. Intra-marrow penetration is useful to ensure some bleeding from the host bone, which provides a supply of osteogenic material to aid in bone regeneration.
- SYRINGE PREPARATION:
 - Pull the piston of the syringe back to be fully extended. Turn the syringe upright (cap upwards), gently tap the container to release the morsels from the cap, and open the cap.
 - Add the patient's blood carefully into the syringe, and close the cap. Do not use marrow in the syringe.
 - Turn the syringe slowly around to wet the morsels completely. If needed, move the piston slowly back and forth to enhance moistening, but be cautious not to pressurize the container excessively. Allow to set until the morsels have soaked up the blood and the mixture has thickened as desired. The mixture can be delivered immediately after the morsels have soaked up the blood, but if the mixture is allowed to set the mixture will become thickened which may be desirable depending on the size and shape of the defect. The maximum recommended immersion time is 60 minutes. Always protect the mixture from contamination by keeping the cap of the syringe on and keeping the syringe inside the sterile field until the mixture has been delivered into the defect.
 - Open the cap before implantation
- PROPER PLACEMENT AND CONTAINMENT OF THE DEVICE:
 - Place the prepared INION BIORESTORE™ bone graft substitute into the intended graft site. The INION BIORESTORE™ should fill the defect as and contact viable bone as much and as tightly as possible. Regeneration will occur best when the material is in tight contact with the bony walls of the defect, and blood and blood vessels can infiltrate the graft material. However, when placing the bone graft substitute into the defect, DO NOT compress the material too strongly to prevent breakage of the material and its porous structure!
 - The INION BIORESTORE™ should be used in areas where the graft can be adequately contained to prevent motion and migration of the material.
 - The INION BIORESTORE™ should fill the bony defect to the highest level of the defect. However, DO NOT overfill the defect! **If the selected implant is too large for the defect, the implant should be shaped to the desired shape according to the above described instructions.**
- SITE CLOSURE:
 - After placement of the INION BIORESTORE™, complete tension-free primary soft tissue closure over the grafted area.
 - DO NOT compromise blood supply to the defect area or blot away the blood/moisture in the placed graft material!
 - DO NOT leave defect open!
 - Degradable or non-degradable membranes or graft containment devices may also be used to achieve closure. The use of a membrane or graft containment device is recommended in cases where the defect is large or limited bony retention is present.
- These instructions are intended as guideline for the use of the INION BIORESTORE™ implants as a part of established techniques. They are not intended to replace or change standard grafting techniques.

Postoperative reminders

- As with any surgical procedure, careful postoperative management is important for optimal healing. Postoperative patient management should follow the same regimen as similar cases utilizing autogenous bone grafting. Standard postoperative practice should be followed.
- Provide the patient with detailed instructions for postoperative care (e.g., regarding appropriate immobilization and hygiene maintenance).

- After treatment of oral and/or dental intraosseous defects, the patient should be cautioned against chewing on or brushing the area for at least 1 to 2 weeks or as long as the clinician deems necessary. Antibacterial rinses should be used during this period.
- Antibiotic and/or analgesic therapy at the discretion of the clinician.
- **INION BIORESTORE™** is MR Safe. The **INION BIORESTORE™** bone graft substitute is entirely made of non-conductive and non-ferrous material components. Measured levels of metallic impurities are negligible.

Evaluation of results

The device is radiopaque: no artifact is created on computer tomography (CT) and magnetic resonance imaging (MRI) by the material.

WARNINGS

- The **INION BIORESTORE™** devices are not intended to withstand the stress of immediate load bearing prior to soft and hard tissue ingrowth without appropriate additional fixation. In cases of fracture fixation or load bearing applications, standard internal or external stabilization techniques must be followed to obtain rigid stabilization in all planes. Instrumentation used in conjunction with the **INION BIORESTORE™** must gain purchase in the host bone. The **INION BIORESTORE™** may not be used to gain screw purchase or to stabilize screw placement. Mechanical loading should not take place for at least 5-6 months after grafting.
- Incorrect selection, placement, or positioning of the implant can cause subsequent undesirable results or breakage of implants. The surgeon should be familiar with bone grafting and internal/external fixation techniques, the devices, the method of application and the surgical procedure prior to performing the surgery.
- Thin soft tissue coverage over the implant or use in areas where the graft cannot be adequately contained (to prevent motion and migration of the material) may increase the risk of complications.

PRECAUTIONS

- DO NOT use for unintended applications! Proper function (i.e., effectivity and safety) of these implants can not be guaranteed in case of off-label use.
- The patient should be warned that the implants can break or loosen as a result of early stress, activity or load bearing to ensure reduced loading to prevent collapse and deformity.

SPECIAL PATIENT POPULATIONS

The **INION BIORESTORE™** has not been clinically tested in pregnant women.

ADVERSE EFFECTS

- Possible complications are the same as to be expected of autogenous bone grafting procedures.
- Premature bending, loosening, breakage or migration of the implants may result from early stress, activity or load bearing.
- Infection can lead to failure of procedure.
- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.
- Complications that may arise as a result of surgery may include: neurovascular injuries due to surgical trauma, superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, delayed union, loss of reduction, failure of fusion, loss of bone graft, graft protrusion and/or dislodgement, tooth sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated tooth root, pain, swelling, and general complications that may arise from anesthesia and/or surgery.

STERILITY

The **INION BIORESTORE™** implants have been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Use only devices that are contained in unopened and undamaged packages. For single use only. Discard remainder. **DO NOT re-sterilize or re-use!** If the device is re-sterilized or re-used, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed. DO NOT use implant beyond the expiration date on the label.

STORAGE

Store at room temperature at a normal relative humidity.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician or dentist.

MANUFACTURER

INION Oy
 Lääkärintäti 2
 33520 Tampere
 Finland
 Tel. +358 10 830 6600
 Fax +358 10 830 6601
 info@inion.com
 www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

Manufacturer

Catalogue number

Lot number

Date of manufacture

Use by

Sterilized using irradiation

Do not re-use

Do not resterilise

Do not use if package is damaged

Keep dry

CE marked product

Read the instructions for use enclosed

Caution, consult accompanying documents

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¹ Präsentation der Society for Biomaterials, 24th Annual Meeting, 1998:511-518

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Inion® on Inion Oy:n rekisteröity tavaramerkki.