ENGLISH

## **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM INSTRUCTIONS FOR USE**

### **DESCRIPTION**

The **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** consists of fixation plates, meshes and screws made of degradable co-polymers composed of L-lactic acid, D-lactic acid, and trimethylene carbonate. These polymers have a long history of safe medical use and they degrade in vivo by hydrolysis into alpha-hydroxy acids that are metabolised by the body. The implants gradually lose their strength during 18-36 weeks in vivo. Bioresorption takes place within two to four years.

The **INION OTPS™** biodegradable implants are offered in different sizes and are designed to be used with customized instrumentation, the **INION®** bone drill bits, bone taps, countersinks, depth gauges, screwdrivers, plate bending pliers and heating devices.

The **INION OTPS™** implants are sterile and non-collagenous.

### **INDICATIONS**

These **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** implants are intended to sustain the relative position of weak bony tissue, e.g., bone grafts, bone graft substitutes, or bone fragments from comminuted fractures. These **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** implants are also indicated for cement restriction in total joint arthroplasty procedures.

Only when used in conjunction with traditional rigid fixation, these **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** implants are intended to sustain the relative position of weak bony tissue in trauma and reconstructive orthopaedic procedures involving the following:

- Long bones
- Flat bones
- Short bones
- Irregular bones
- Appendicular skeleton
- Thorax

When used alone (without traditional rigid fixation), these **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** implants are intended to maintain the relative position of bone grafts or bone graft substitutes in reconstructive orthopaedic procedures involving:

- Tumor resections where bone strength has not been compromised.
- Iliac crest harvest sites.

### **CONTRAINDICATIONS**

These **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** implants are not intended for use in and are contraindicated for:

- Spine indications.
- Load bearing indications unless used in conjunction with traditional rigid fixation.
- Active or potential infection.
- Patient conditions including limited blood supply, insufficient quantity or quality of bone; and where patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse).

### **INFORMATION FOR USE**

#### **Surgical considerations and reminders**

- Prophylactic perioperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.
- Maintain sterile field throughout the procedure.
- Proper exposure using standard surgical procedure.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Good alignment/reduction of the fracture/osteotomy.
- Radiographs can be taken to check the alignment/reduction after fixation.
- Meticulous hemostasis and complete primary soft tissue closure over the implant are essential.

#### **Preparation of the meshes**

The meshes should be contoured to the bone according to these instructions (see also the instructions for use of the **INION THERMO+™** water bath):

- Place a sterile **INION THERMO™** drape to the **INION THERMO+™** water bath. DO NOT switch on the water bath until filled with liquid!
- Pour sterile water or saline (min 0.4 litres max 0.7 litres) into the water bath.
- Switch the water bath on. Wait until the liquid heats up to the setting temperature 70°C / 158°F, this takes about 30 minutes.
- Choose appropriate **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** mesh for the indication.
- Immerse the mesh in the water bath for approximately one minute. After water bath treatment: The mesh is most malleable for 10-15 seconds, can be shaped as desired and may be cut to the desired size with standard surgical scissors.
- Maximum total immersion time of the plate/mesh in the water bath is 30 minutes. Prolonged immersion may weaken the plate/mesh.
- Place the mesh against bone at its intended location. Contour the mesh to desired shape using fingers and/or the **INION®** plate bending pliers.
- If contouring is not satisfactory, the water bath treatment can be repeated. If only part of the mesh needs further adaptation dip only that portion of the mesh or bend once cooled.
- Repeated bending can weaken the mesh, especially once cooled. Warm bending can be carried out up to three times.
- Take care not to perforate the **INION THERMO™** drape with surgical instruments.
- Take care to ensure that the mesh is accurately adapted to the bone contours.
- DO NOT use mesh(es) that show signs of damage!

- DO NOT activate the screws by heating in the water bath!

### Fastening of the meshes

- Choose appropriate **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** screws for fixation of the mesh(es). If needed, the screws can be cut to desired length with standard surgical scissors.

**NOTE: Always rinse the drill hole before screw insertion.**

#### • DRILL & MANUAL TAP METHOD:

- Using the appropriate (corresponding to the screw diameter) **INION®** bone drill bit attached to a slow speed drill (maximum speed of 2000 rpm), drill a suitable hole for the screw into the bone. Use irrigation. Always drill the holes for the screws through the holes of the mesh to ensure correct location and mesh-hole size. Take care to hold the bone drill parallel to the long axis of the drill hole while drilling to prevent drill-wobble induced enlargement of the hole.
- Tap the screw holes through the holes of the mesh manually using the appropriate (corresponding to the screw diameter) **INION®** bone tap before screw insertion. Hold the bone tap parallel to the long axis of the drill hole both during insertion and reversing of the bone tap. Only use gentle pressure when tapping to preserve the prepared threads. Care must be taken to tap to a sufficient depth to allow complete seating of the screw. The prepared hole should be irrigated prior to screw insertion to flush out bone debris.

#### • SCREW LOADING & PLACEMENT:

- Align the appropriate (corresponding to the screw diameter) **INION®** screwdriver with the top of the screw. Turn the screwdriver until the screw rotates. Gently push the screwdriver into the screw head and lift the screw out rotating it slightly. Take care to ensure that the screwdriver is aligned properly before pushing it into the screw, since there is a risk of damaging the plastic screw head.
- Hold the screwdriver and the screw parallel to the long axis of the drill hole and while keeping the mesh in position, insert the screw fully into the drill hole. DO NOT overtighten!

### Postoperative reminders

- As with any surgical procedure, careful postoperative management is important for optimal healing.
- Provide the patient with detailed instructions for postoperative care (e.g., regarding immobilization and hygiene maintenance).
- Use appropriate additional immobilization (e.g., a suitable cast, brace and/or crutches) during bone healing.
- Antibiotic therapy at the discretion of the clinician.

### Evaluation of results

Radiographs can be taken to evaluate healing.

### WARNINGS

- The **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** implants provide fixation and are not intended to replace healthy bone or withstand the stress of full load bearing.
- Incorrect selection, placement, positioning, and fixation of the implant can cause subsequent undesirable results. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.

### PRECAUTIONS

- Instruments are available to aid accurate implantation of the **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM**. Surgical instruments are subject to wear with normal usage and may break. Surgical instruments are only to be used for their intended purpose. All instruments are to be regularly inspected for wear and damage. Use only the **INION®** instruments.
- Use only **INION OTPS™** screws to secure **INION OTPS™** meshes and plates. DO NOT use **INION OTPS™** screws to secure plates and meshes other than those of **INION®** (e.g., metal plates).
- DO NOT use for unintended applications! Proper function (i.e., effectivity and safety) of these implants can not be guaranteed in case of off-label use.
- The patient should be warned that the implants can break or loosen as a result of early stress, activity or load bearing. Premature discontinuation of the cast or other immobilization technique may cause non-union or mal-union of the fracture or osteotomy.

### ADVERSE EFFECTS

Complications are similar to those encountered with any other method of internal fixation:

- Premature bending, loosening, fracture or migration of the devices may result from early stress, activity or load bearing.
- Infection can lead to failure of the procedure.
- Neurovascular injuries can occur due to surgical trauma.
- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.

### STERILITY

The **INION OTPS™ BIODEGRADABLE FIXATION SYSTEM** implants have been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Use only devices that are contained in unopened and undamaged packages. For single use only. DO NOT re-sterilize. If the device is re-sterilized, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed. DO NOT use implant beyond the expiration date on the label.

### STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F) at a normal relative humidity. Product should not exceed a maximum temperature defined on the indicator. Do not use if the indicator dot is black.

### CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician.

### MANUFACTURER


INION Oy

Lääkärintie 2  
33520 Tampere  
Finland  
Tel. +358 10 830 6600  
Fax +358 10 830 6601  
info@inion.com  
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

## SYMBOLS USED IN LABELING

	Manufacturer		
	Catalogue number		Lot number
	Date of manufacture		Use by
	Sterilized using irradiation		CE marked product
	Storage temperature range		Keep dry
	Caution, consult accompanying documents		Do not re-use

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Inion® on Inion Oy:n rekisteröity tavaramerkki.