ENGLISH

INION OTPS™ BIODEGRADABLE PIN INSTRUCTIONS FOR USE

DESCRIPTION

The **INION OTPS™ BIODEGRADABLE PIN** is intended to maintain accurate alignment of fragments of fractured bone in the presence of appropriate immobilization.

The **INION OTPS™ BIODEGRADABLE PIN** is made of degradable co-polymers composed of L-lactic acid, D-lactic acid and trimethylene carbonate. These polymers have a long history of safe medical use and they degrade in vivo by hydrolysis into alpha-hydroxy acids that are metabolised by the body. The coloured pins are dyed green for better visualization during the surgical procedure by a minimal amount of Drug and Cosmetic (D&C) Green No. 6, which is used in several biodegradable sutures. The pin gradually loses its strength during 18-36 weeks. Bioresorption takes place within two to four years.

The **INION OTPS™ BIODEGRADABLE PIN** is offered in different sizes and it is designed to be used with customized instrumentation, e.g. the **INION OTPS™ TROCAR TIP DRILL**, **PIN APPLICATORS**, and **PERCUTANEOUS** and **ARTHROSCOPIC** instruments.

The **INION OTPS™ BIODEGRADABLE PIN** is sterile and non-collagenous.

INDICATIONS

The **INION OTPS™ BIODEGRADABLE PIN** is indicated for maintenance of alignment and fixation of bone fractures, osteotomies, arthrodeses or bone grafts in the presence of appropriate additional immobilization (e.g. rigid fixation implants, cast, brace).

CONTRAINDICATIONS

The **INION OTPS™ BIODEGRADABLE PIN** should not be used in fractures and osteotomies of diaphyseal bone or in cases with insufficient quality or quantity of bone. Other contraindications are active or potential infections, patient conditions including limited blood supply, and where patient cooperation cannot be guaranteed (e.g. alcoholism, drug abuse). There are currently no known additional contraindications to the use of the **INION OTPS™ BIODEGRADABLE PIN**.

INFORMATION FOR USE

Surgical considerations and reminders

- Prophylactic perioperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.
- Maintain sterile field throughout the procedure.
- Proper exposure using standard surgical procedure.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Good alignment/reduction of the fracture/osteotomy followed by fixation with clamp(s).
- Radiographs can be taken before wound closure to check the alignment/reduction after fixation.
- Meticulous hemostasis and complete primary skin closure over the implant are essential.

Implantation of the pin with reusable instruments

- Choose appropriate **INION OTPS™ BIODEGRADABLE PIN** for the indication.
- Using the appropriate (corresponding to the pin diameter) **INION OTPS™ TROCAR TIP DRILL**, drill a hole through the fracture/osteotomy plane.
- Irrigate the drill hole.
- Introduce the pin into the hole by hand or using the appropriate (corresponding to the pin diameter) **INION OTPS™ PIN APPLICATOR**. The pin should sink slightly into the drilled hole. If not, widen the opening of the hole a little with the drill.
- Push the pin applicator onto the pin and press the applicator shaft to contact with the pin. The **INION®** instruments must be used to ensure proper insertion of the pin into the drill hole. During insertion of the pin, hold the applicator and the pin parallel to the long axis of the drill hole to prevent bending of the pin.
- Tap the applicator shaft into the cylinder with a small mallet so that the entire pin is forced fully into the drill hole. This prevents the head of the pin protruding which could cause soft tissue irritation.
- When necessary, a hot wire or scissors can be used to cut the pin before insertion or after insertion if the pin is too long.
- Two or more drill holes with fixing pins can be applied if necessary (depending on the nature and size of the fracture).
- In case a pin breaks during insertion, the broken pin can be drilled out and a new pin can be placed into the same drill hole.
- If the end result of fracture reduction or fixation is not satisfactory, the pin(s) and fixation can be removed by drilling the pin(s) out.
- When performing percutaneous pin fixation, use the **INION OTPS™ PERCUTANEOUS** instruments.
- When performing arthroscopic pin fixation, use the **INION OTPS™ ARTHROSCOPIC** instruments.

Implantation of the pin with disposable single-use instruments included in the INION OTPS™ Pin Kit

- Using the K-wire provided in the **INION OTPS™ Pin Kit**, drill a hole through the fracture/osteotomy plane.
- Irrigate the drill hole.
- Introduce the hooked end of the depth gauge provided in the kit into the entire depth of the drill hole. If the drill hole passes entirely through the bone, pull back on the depth gauge until the hook catches the opposite cortex. Slide the indicator down the depth gauge until it contacts the surface of the bone. Thereafter remove the depth gauge from the drill hole and take care not to change the position of the stop.
- Use scissors or a hot wire to cut the pin to match the distance between the hooked end and the stop of the depth gauge.
- Introduce the pin into the cannula provided in the kit.
- Partially insert the introducer provided in the kit into the cannula so the tip of the pin protrudes from the cannula. Introduce the pin tip into the drilled hole. The pin should sink slightly into the drill hole. If not, widen the opening of the hole a little with the K-wire.

- Alternatively, the pin can be introduced into the hole by hand and then the cannula placed over the pin, following which the introducer is placed into the cannula
- Press the introducer to contact with the pin. During insertion of the pin, hold the cannula and the pin parallel to the long axis of the drill hole to prevent bending of the pin.
- Tap the introducer into the cannula with a small mallet so that the entire pin is forced fully into the drill hole. This prevents the head of the pin protruding which could cause soft tissue irritation.
- When necessary, a hot wire or scissors can be used to cut the pin after insertion if the pin is too long.
- Two or more drill holes with fixing pins can be applied if necessary (depending on the nature and size of the fracture).
- In case a pin breaks during insertion, the broken pin can be drilled out and a new pin can be placed into the same drill hole.
- If the end result of fracture reduction or fixation is not satisfactory, the pin(s) and fixation can be removed by drilling the pin(s) out.
- When performing percutaneous pin fixation, use the **INION OTPS™ PERCUTANEOUS** instruments.
- When performing arthroscopic pin fixation, use the **INION OTPS™ ARTHROSCOPIC** instruments.

Post operative reminders

- As with any surgical procedure, careful postoperative management is important for optimal healing.
- Provide the patient with detailed instructions for postoperative care.
- Use appropriate additional immobilization (e.g. a suitable cast, brace and/or crutches) during bone healing.
- Antibiotic therapy at the discretion of the clinician.

Evaluation of results

Radiographs can be taken to evaluate bone healing.

WARNINGS

- The **INION OTPS™ BIODEGRADABLE PIN** implants provide fixation and are not intended to replace healthy bone or withstand the stress of full load bearing.
- Incorrect selection, placement, positioning, and fixation of the implant can cause subsequent undesirable results. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.

PRECAUTIONS

- Instruments are available to aid accurate implantation of the **INION OTPS™ BIODEGRADABLE PIN**. Surgical instruments are subject to wear with normal usage and may break. Surgical instruments are only to be used for their intended purpose. All reusable instruments are to be regularly inspected for wear and damage. All single-use instruments should be discarded after use! Use only the **INION OTPS™ PIN** instruments.
- When two or more pins are used, insert the pins at divergent angles to one another rather than parallel, for best results. When using several pins, firstly insert the corresponding K-wires (Trocar Tip Drills) through the fragment and both cortices to secure the fragment. Then remove one of the K-wires, irrigate and measure the depth of the drill hole, and insert an appropriate length pin of the same diameter as the K-wire used into the drill hole. Thereafter, replace the other K-wires one by one in the same manner.
- DO NOT use for unintended applications! Proper function (i.e., effectivity and safety) of these implants can not be guaranteed in case of off-label use.
- The patient should be warned that premature bending, loosening, breakage or migration of the pin may result from early weight bearing, stress and activity.

SPECIAL PATIENT POPULATIONS

The effect of the **INION OTPS™ BIODEGRADABLE PIN** upon the healing of growth plate has not been tested clinically.

ADVERSE EFFECTS

Complications are similar to those with any method of internal fixation:

- Premature bending, loosening, breakage or migration of the devices may result from early stress, activity or load bearing.
- Infection can lead to failure of the procedure.
- Neurovascular injuries can occur due to surgical trauma.
- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.

STERILITY

The **INION OTPS™ BIODEGRADABLE PIN** implants have been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Use only devices that are contained in unopened and undamaged packages. For single use only. DO NOT re-sterilize. If the device is re-sterilized, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed. DO NOT use implant beyond the expiration date on the label.

STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F) at a normal relative humidity. Product should not exceed a maximum temperature defined on the indicator. Do not use if the indicator dot is black.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician.

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
Fax +358 10 830 6601
info@inion.com
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

	Manufacturer		
	Catalogue number		Lot number
	Date of manufacture		Use by
	Sterilized using irradiation		CE marked product
	Storage temperature range		Keep dry
	Caution, consult accompanying documents		Do not re-use

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INION Oy
Inion® on Inion Oy:n rekisteröity tavaramerkki.