12000431 (5/2019-04-01)

ENGLISH

INION OTPS FREEDOMPLATE™ INSTRUCTIONS FOR USE

DESCRIPTION

The **INION OTPS FREEDOMPLATE™** biodegradable implants are made of degradable co-polymers composed of L-lactic acid, D-lactic acid and trimethylene carbonate. These polymers have a long history of safe medical use and they degrade *in vivo* by hydrolysis into alpha-hydroxy acids that are metabolised by the body. The implants gradually lose their strength during 18-36 weeks. Bioresorption takes place within two to four years.

The **INION OTPS FREEDOMPLATE™** implants are offered in different sizes and are designed to be used with customized instrumentation, e.g., the **INION®** bone drill bits, bone taps, drill taps, countersinks, drill guides, depth gauges, screw cutting pliers, screw driver, plate bending pliers, temperature cautery instruments and heating devices.

The **INION OTPS FREEDOMPLATE™** implants are sterile and non-collagenous.

INDICATIONS

A. General indications: In the presence of appropriate additional immobilization or fixation, indicated for maintenance of alignment and fixation of bone fractures, osteotomies, arthrodeses or bone grafts, and maintenance of relative position of weak bony tissue (e.g., bone grafts, bone graft substitutes, or bone fragments from comminuted fractures), in trauma and reconstructive procedures.

B. Specific indications: Long bones, flat bones, short bones, irregular bones, appendicular skeleton, and thorax.

CONTRAINDICATIONS

The **INION OTPS FREEDOMPLATE™** implants are not intended for use in and are contraindicated for:

- Active or potential infection.
- Patient conditions including limited blood supply, insufficient quantity or quality of bone; and where patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse).
- Load bearing applications unless used in the presence of appropriate additional immobilization or fixation.
- Spinal indications.
- Mandibular tumor resection.
- Mandible unless used with appropriate maxillomandibular fixation.

INFORMATION FOR USE

Surgical considerations and reminders

- The operating principles follow those used for other methods of internal fixation as described in textbooks.
- Prophylactic preoperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.
- Maintain sterile field throughout the procedure.
- Proper exposure using standard surgical procedure.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Good alignment/reduction of the fracture/osteotomy.
- Radiographs can be taken before wound closure to check the alignment/reduction after fixation.
- Meticulous hemostasis and complete primary soft tissue closure over the implant are essential.

Preparation of the plates

The plates should be contoured to the bone according to these instructions (see also the instructions for use of the **INION THERMO+™** water bath):

- Place a sterile **INION THERMO™** drape to the **INION THERMO+™** water bath. DO NOT switch on the water bath until filled with liquid!
- Pour sterile water or saline (min 0.4 litres max 0.7 litres) into the waterbath.
- Switch the water bath on. Wait until the liquid heats up to the setting temperature 70°C / 158°F, this takes about 30 minutes.
- Choose appropriate **INION OTPS FREEDOMPLATE™** for the indication.
- Immerse the plate in the water bath for approximately one minute. After water bath treatment: The plate is most malleable for 10-15 seconds, can be shaped as desired and may be cut to the desired size with standard surgical scissors or cutting pliers (easiest immediately after picking up the plate from the water bath).
- Maximum total immersion time of the plate/mesh in the waterbath is 30 minutes. Prolonged immersion may weaken the plate/mesh.
- **Take care to cut the plate to appropriate dimensions corresponding to the implantation site and plate size typical for the intended use.** DO NOT implant too small or unnecessary large plates for the purpose! When cutting the plate, DO NOT cut too close to the pilot holes where you are planning to insert screws. At minimum, half of the distance between the original pilot holes should remain around any pilot hole intended for screw placement. Take care to round any sharp edges or corners of the trimmed plate before implantation to prevent irritation of the soft tissue around the plate.
- Place the water bath treated plate against bone at its intended location.
- When using concave **INION OTPS FREEDOMPLATE™**, ensure that the concave side (not the convex side) of the plate is placed against bone.
- Contour the plate to desired shape by hand and/or by using the **INION®** plate bending pliers.
- If contouring or shape of the plate is not satisfactory, the water bath treatment can be repeated by immersing the plate in the water bath for additional 20-30 seconds. If only part of the plate needs further adaptation dip only that portion of the plate or bend once cooled.

- Repeated bending can weaken the plate, especially once cooled. Warm bending can be carried out up to three times. Please note that the plate returns to its original shape in the water bath.
- Take care not to perforate the **INION THERMO™** drape with surgical instruments.
- Take care to ensure that the plate is accurately adapted to the bone contours.
- DO NOT use implants that show signs of damage!
- DO NOT activate the screws by heating in the water bath!

Fastening of the plates

- Choose appropriate **INION CPS®/OTPS™ Ø2.5/2.8/3.1 mm** screws or **INION FreedomScrew™ Ø2.7/3.5/4.0/4.5 mm** products for fixation of the plate. If needed, the screws can be cut to desired length with the **INION®** screw cutting pliers. DO NOT use **INION CPS®/OTPS™ Ø2.0 mm** screws or **INION FreedomScrew™ Ø2.0 mm** products as a primary fixation screw for the **INION OTPS FREEDOMPLATE™**.
- When using **INION FreedomScrew™ Ø2.7/3.5/4.0/4.5 mm** products for fastening the **INION OTPS FREEDOMPLATE™**, please refer also to the product line-specific instructions for use.
- Always insert screws through the pilot holes of the plate closest to the fracture/osteotomy line. For optimal fixation strength the screws should be inserted bicortically. To ensure that plate edges remain in tight contact with bone, place screws close to the edges of the plate. To achieve an optimal screw placement, the plate is cut exactly in between the pilot holes. **NOTE: Always rinse the drill hole before screw insertion. If needed, use drill bit, bone taps and drill taps through appropriate drill guide. This is to ensure correct hole position and direction, to evaluate appropriate tapping depth and to protect the surrounding tissues from the rotating instrument.**
- Cannulated instruments must be used with cannulated screws. It is recommended to make an x-ray check of the location of K-wire before drilling.
- **DRILL & MANUAL TAP METHOD:**
 - The material of the plate allows drilling of holes for screws through the pilot holes of the plate. Please note that drilling of holes between the pilot holes may weaken the mechanical properties of the plate. Always use irrigation when drilling through the plate (important in order to avoid plate hole enlargement due to melting the plate material! Enlargement of the plate hole may compromise screw plate stability and interlocking strength!).
 - Place the tip of an appropriate (corresponding to the screw diameter) **INION®** drill guide by the intended pilot hole of the plate, angulate the drill guide in the desired direction for the screw, and maintain the position and direction of the drill guide thereafter. However, take care not to angulate the drill guide too much to avoid drilling through the other screws or into a joint. Angulation of the screw may also cause protrusion of the screw head.
 - If needed, an appropriate **INION®** countersink instrument can be used to enable more low-profile seating of the screw head or the head of the screw can simply be cut off along the plate surface with cutting pliers or 1-2 mm above the plate surface with a temperature cautery to prevent secondary soft tissue irritation.
 - Using the appropriate **INION®** drill bit attached to a slow speed drill (maximum speed of 2000 rpm), drill a suitable hole for the screw through the intended pilot hole of the plate into the bone. This is important in order to avoid plate hole enlargement due to melting the plate material! Enlargement of the plate hole may compromise screw plate stability and interlocking strength. Use irrigation. Always drill the holes for the screws through the drill guide and the plate to ensure correct final plate-hole size, drill-hole position and direction and to protect the surrounding tissues from the drill. Avoid drilling too deep or into joints to prevent surgical trauma. Hold the drill guide and the bone drill parallel to the long axis of the drill hole while drilling to prevent drill-wobble induced enlargement of the bone or plate hole. Enlargement of the bone or plate hole may compromise screw plate stability and interlocking strength, and bone fixation strength!. During drilling, press the plate tightly against bone.
 - Always measure the depth of the drill hole with a depth gauge through the hole of the plate and choose screw length accordingly. If needed, the screws can be cut to desired length with the **INION®** screw cutting pliers. If the intention is to cut the head of the screw off after screw insertion, the screw length can be selected to be longer than the depth of the drill hole.
 - Before screw insertion tap the drill hole manually by using the appropriate (corresponding to the screw diameter) **INION®** bone tap, and always insert the bone tap into the drill hole through the corresponding hole of the plate. If needed, use bone taps through an appropriate drill guide. This is to ensure correct hole position and direction, to evaluate appropriate tapping depth, and to protect the surrounding tissues from the rotating instrument. During bone tap insertion, press the plate tightly against bone. Hold the bone tap parallel to the long axis of the drill hole both during insertion and reversing of the bone tap. This is important in order to avoid bone and plate hole enlargement! Enlargement of the bone or plate hole may compromise screw plate stability and interlocking strength, and bone fixation strength. Use irrigation. Only use gentle pressure when tapping to preserve the prepared threads. Care must be taken to tap to a sufficient depth to allow complete seating of the screw. **NOTE: Always rinse the drill hole before screw insertion to flush out bone debris and to enable complication-free screw insertion!**
 - When drilling and tapping is carried out according to the instructions above, the plate holes will also become threaded and the screws will interlock with the plate. This interlocking provides sufficient screw plate stability and fixation strength enabling the head of the screw to be cut off at the plate surface with cutting pliers or 1-2 mm above the plate surface with a temperature cautery after screw insertion to reduce the risk of soft tissue irritation due to protruding screw head. Take care to round any sharp edges or corners of the cut screw head to prevent irritation of the soft tissue over the plate.
 - Please note that if the head of the screw is not cut off, the appropriate **INION®** countersink instrument can be used prior to screw insertion in order to make space for the screw head for more low-profile seating of the screw head to avoid secondary soft tissue irritation (this should be done especially if the plate is implanted under a thin soft tissue layer).
 - Please note that before screw insertion, some of the pilot holes of the plate can be drilled larger by using a drill bit with equal outer diameter with the selected screw to create unthreaded non-interlocking holes that allow compression of the plate against the bone surface by tightening the screws. Alternatively the pilot hole can be enlarged with appropriate **INION®** countersink instrument. In any case, take care to enlarge only the hole of the plate, not the hole in the bone!
Please note that the head of the screw must not be cut off if the hole of the plate is made non-interlocking by enlarging the hole as described above.
- **DRILL TAP METHOD (with INION FreedomScrew™ products):**
 - Instruments for simultaneous drilling and tapping are marked with a yellow stripe on the instrument shaft.
 - The material of the plate allows creating threaded holes for screws through the pilot holes of the plate. Please note that creating holes between the pilot holes may weaken the mechanical properties of the plate. **NOTE: Always use irrigation when drilling/tapping through the plate.** This is important in order to avoid plate hole enlargement due to melting the plate material! Enlargement of the plate hole may compromise screw plate stability and interlocking strength.

- Place the tip of an appropriate (corresponding to the screw diameter) **INION®** drill guide by the intended pilot hole of the plate. Angulate the drill guide in the desired direction for the screw, and maintain the position and direction of the drill guide thereafter. However, take care not to angulate the drill guide too much to avoid drilling/tapping through the other screws or into a joint. Angulation of the screw may also cause protrusion of the screw head.
- Use an appropriate **INION®** drill tap instrument attached to a slow speed drill, which has a reverse function, to create a suitable, threaded hole for the selected **INION FreedomScrew™** through the intended pilot hole of the plate into the bone. This is important in order to avoid plate hole enlargement due to melting the plate material! Enlargement of the plate hole may compromise screw plate stability and interlocking strength. Use irrigation. Advance the appropriate **INION®** drill tap slowly using a moderate amount of pressure to push the instrument forward against the plate and bone in the direction parallel to the long axis of the bone hole being created. Always create the holes for the screws through the drill guide and the plate. This ensures the correct size of the final plate and the hole, and the correct position and direction of the hole. This also helps to evaluate appropriate tapping depth, and to protect the surrounding tissues from the rotating instrument. Avoid advancing the drill tap instrument too deep or into joints to prevent surgical trauma. Take care to hold the drill guide and the drill tap instrument parallel to the long axis of the bone hole when inserting and when reversing the drill tap instrument during plate and bone hole creation. This prevents plate or bone hole enlargement induced by instrument wobble, or inappropriate threads on plate or bone. Enlargement of the bone or plate hole or inappropriate threads may compromise screw plate stability and interlocking strength, as well as bone fixation strength.

During the insertion of the drill tap instrument and the creation of the plate and bone hole, press the plate tightly against the bone. To remove the drill tap from the bone, set the drill machine to reverse first. To preserve the prepared threads, pull gently when reversing the drill tap instrument from the threaded bone hole by using a slow speed drill. Carefully create a sufficient depth for bone hole to allow complete seating of the screw. DO NOT advance drill tap instrument deeper into the bone than its threaded part allows in order to preserve the prepared threads! **NOTE: Always rinse the prepared hole before screw insertion to flush out bone debris and to enable complication-free screw insertion!**

- If needed, the screws can be cut to desired length with the **INION®** screw cutting pliers. If the intention is to cut the head of the screw off after screw insertion, you can select a screw that is longer than the depth of the bone hole.
- When the above described DRILL TAP METHOD is used in the fixation of plates, the recommendation is to always cut off the screw head 1-2 mm above the plate surface with a temperature cautery after screw insertion. This ensures optimal fixation strength and removes or rounds any sharp edges of the plate drill hole or cut screw head. Sharp edges can cause soft tissue irritation.
- If the screw head is not cut off, the appropriate **INION®** countersink instrument can be used prior to screw insertion in order to make space for a more low-profile seating. This is to avoid secondary soft tissue irritation. When using a countersink instrument, take care to enlarge only the hole of the plate, not the hole in the bone! **NOTE: The head of the screw must not be cut off if the hole of the plate has been enlarged by the countersink instrument, and thus the plate hole is made non-interlocking.**

• SCREW LOADING & PLACEMENT:

- Attach the appropriate **INION®** or other compatible screwdriver shaft with the appropriate holding sleeve. Align the sleeve with an insertion adapter of the screw and connect the sleeve with the insertion adapter by sliding the sleeve towards the screw. Take care to ensure that the sleeve is aligned and connected properly with insertion adapter.
- Hold the screw driver and the screw parallel to the long axis of the drill hole and while keeping the plate in position tightly against bone surface, insert the screw through the hole of the plate fully into the drill hole. DO NOT overtighten!
- If the countersink instrument has not been used for low-profile seatings for the screw heads, after successful screw insertions, cut off the screw heads. This prevents soft tissue irritation as described in sections DRILL & MANUAL TAP METHOD and DRILL TAP METHOD above.
- If a screw gets stuck during its insertion and cannot be fully inserted, remove the screw from the drill hole by reversing the screw out by using the screw driver. If the screw breaks or the head of the screw gets damaged during screw insertion or reversing and the screw cannot be removed with the screw driver, drill through the screw with the drill bit that was originally used to create the drill hole for the screw. Alternatively the protruding part of the screw can simply be cut off if it is estimated that the screw is deep enough for adequate fixation. Always use the drill guide when drilling through screws and take care to angulate the drill guide and the drill parallel to the long axis of the screw (Enlargement of the bone or plate hole may compromise screw plate stability and interlocking strength, and bone fixation strength!).
- After screw removal by reversing the screw with the screw driver:
 1. Rinse the drill hole and measure the length of the drill hole again to check the correct screw length. If the previously selected screw length was longer than the length of the bone tunnel, select a shorter screw (corresponding to the length of the drill hole) or cut the previously selected screw (if undamaged) to correct length with the **INION®** screw cutting pliers, or re-drill the hole if a longer screw hole length or bicortical fixation is desired. If the originally selected screw length was equal to or shorter than the length of the drill hole, re-tap the drill hole with the previously used **INION®** bone tap (corresponding to the screw diameter) following the above described instructions. Care must be taken to insert the bone tap deep enough to create threads along the whole length of the bone tunnel and fully through all involved cortices to allow complete seating of the screw. Always use the drill guide when re-drilling, and take care to angulate the drill guide and the instruments parallel to the long axis of the original drill hole. Re-drilling or re-tapping in a different angle than that chosen for the original hole may damage the previously created threads and cause enlargement of the hole or it may cause plate hole enlargement! Enlargement of the bone or plate hole may compromise screw plate stability and interlocking strength, and bone fixation strength!
 2. Rinse the drill hole before screw insertion to flush out bone debris and to enable complication-free screw insertion.
 3. Insert the appropriate **INION CPS®/OTPS™** screw or **INION FreedomScrew™** product following the above described instructions for screw insertion. The previously removed screw can be used if it is undamaged.
- After screw removal by drilling through the biodegradable screw:
 1. Rinse the drill hole, measure the length of the hole, and re-tap the hole with the previously used **INION®** bone tap (corresponding to the screw diameter).
 2. Rinse the drill hole before screw insertion to flush out bone debris and to enable complication-free screw insertion.
 3. Insert a new appropriate **INION CPS®/OTPS™** screw or **INION FreedomScrew™** product (corresponding to the measured drill hole length and original screw diameter) into the drill hole following the above described instructions.
- If a screw feels loose, remove the screw from the drill hole by reversing the screw out. Use the screwdriver or drill through the screw with the drill bit that was originally used for creating the drill hole for the screw. Re-tap the drill hole with one size larger (corresponding to the screw diameter) **INION®** bone tap or drill tap, and insert a corresponding larger diameter **INION CPS®/OTPS™** screw or **INION FreedomScrew™** product into the drill hole following the above described instructions.

Post operative reminders

- As with any surgical procedure, careful postoperative management is important for optimal healing.
- Provide the patient with detailed instructions for postoperative care (e.g., regarding appropriate immobilization, hygiene maintenance, and/or soft non chewing diet if needed).
- Use appropriate additional immobilization (e.g., a suitable plaster cast; brace; orthosis; crutches; or maxillomandibular fixation, especially when treating mandibular fractures and bimaxillary osteotomies, e.g., guiding elastics) or fixation during bone healing.
- Antibiotic therapy at the discretion of the clinician.

Evaluation of results

Radiographs, magnetic resonance imaging (MRI) and/or computer tomography (CT) scans can be taken postoperatively to evaluate healing. The device is MR Safe.

WARNINGS

- The **INION OTPS FREEDOMPLATE™** implants provide fixation and are not intended to replace healthy tissues or withstand the stress of full load bearing.
- Incorrect selection, placement, positioning, or fixation of the implant can cause subsequent undesirable results or breakage of implants or instruments. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.
- The implants can break or loosen if the plate is used as a lever to reduce a fracture/osteotomy!
- Thin soft tissue coverage over the implant may increase the risk of complications.
- Implant may be covered with fibrous tissue capsule when implanted.

PRECAUTIONS

- Instruments are available to aid accurate implantation of the **INION OTPS FREEDOMPLATE™** implants. Surgical instruments are subject to wear with normal usage and may break. Surgical instruments are only to be used for their intended purpose. All instruments are to be regularly inspected for wear and damage. • Use only the **INION CPS®/OTPS™ Ø2.5/2.8/3.1 mm** screws or **INION FreedomScrew™ Ø2.7/3.5/4.0/4.5 mm** products to secure the **INION OTPS FREEDOMPLATE™**. DO NOT use the **INION CPS®/OTPS™** screws or **INION FreedomScrew™** products to secure plates or sheets other than those of **INION®** (e.g., metal plates or sheets).
- DO NOT use for unintended applications! Proper function (i.e., effectivity and safety) of these implants can not be guaranteed in case of off-label use.
- The patient should be warned that the implants can break or loosen as a result of early stress, activity or load bearing. Premature discontinuation of additional postoperative immobilization or fixation may cause non-union or mal-union.

SPECIAL PATIENT POPULATIONS

Patients with diabetes, bruxism, and temporomandibular joint malfunction or malformation may have an increased risk of complications following fixation in the maxilla and mandible. Patients with symptoms of sleep apnea have increased risk for complications following fixation of osteotomies in the maxilla and mandible.

ADVERSE EFFECTS

Complications are similar to those with any method of internal fixation:

- Premature bending, loosening, breakage or migration of the devices may result from early stress, activity or load bearing.
- Infection can lead to failure of procedure.
- Neurovascular injuries can occur due to surgical trauma.
- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.

STERILITY

The **INION OTPS FREEDOMPLATE™** implants have been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Use only devices that are contained in unopened and undamaged packages. For single use only. DO NOT resterilize or reuse! If the device is re-sterilized or reused, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed. DO NOT use implant beyond the expiration date on the label.

STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F) at a normal relative humidity. Product should not exceed a maximum temperature defined on the indicator. Do not use if the indicator dot is black.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician.

MANUFACTURER

INION Oy
Lääkärintäti 2
33520 Tampere
Finland
Tel. +358 10 830 6600
Fax +358 10 830 6601
info@inion.com
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

Manufacturer

Catalogue number

Date of manufacture

Sterilized using irradiation

Storage temperature range

Caution, consult accompanying documents

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Inion® on Inion Oy:n rekisteröity tavaramerkki.