ENGLISH

## INION FreedomScrew™ INSTRUCTIONS FOR USE

### DESCRIPTION

The **INION FreedomScrew™** products are fixation screws made of degradable co-polymers composed of L-lactic acid and D-lactic acid. These polymers have a long history of safe medical use and they degrade in vivo by hydrolysis into alpha-hydroxy acids that are metabolised by the body. Based on in vitro testing, the screws retain their initial strength up to 12 weeks after implantation and gradually lose their strength thereafter. Bioresorption takes place within two to four years.

The **INION FreedomScrew™** products are offered in different sizes and are primarily designed to be used with customized instrumentation, the **INION®** bone drill bits, bone taps, drill taps, holding sleeves, screwdriver shafts, drill sleeves, countersinks, depth gauges and temperature cautery instruments.

The **INION FreedomScrew™** products are also designed to be used with ISO 9714 compatible instrumentation. Corresponding **INION FreedomScrew™** and instrument sizes are presented in the Table 1 (below).

The **INION FreedomScrew™** products are delivered with the disposable single-use metallic insertion adapter with a standard hex socket.

The **INION FreedomScrew™** products are sterile and non-collagenous.

### INDICATIONS

The **INION FreedomScrew™** products are intended for maintenance of alignment and fixation of bone fractures, comminuted fractures, osteotomies, arthrodeses or bone grafts (i.e., autografts or allografts) in the presence of appropriate additional immobilization (e.g., rigid fixation implants, cast or brace).

In addition, the **INION FreedomScrew™ Ø3.5/4.0/4.5 mm** products are specifically intended for use in following indications:

- A. General indications: maintenance of reduction and fixation of cancellous bone fractures, osteotomies or arthrodeses of the upper extremity, ankle and foot in the presence of appropriate brace and/or immobilization.
- B. Specific indications: fractures and osteotomies of the malleoli, and ankle fractures.

### CONTRAINDICATIONS

The **INION FreedomScrew™** products are not intended for use in and are contraindicated for:

- Active or potential infection
- Patient conditions including limited blood supply, insufficient quantity or quality of bone; and where patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse)
- Spine indications
- High-load bearing indications (e.g., diaphyseal fractures of long bones) unless used in conjunction with traditional rigid fixation

### INFORMATION FOR USE

#### Surgical considerations and reminders

- Prophylactic perioperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.
- Maintain sterile field throughout the procedure.
- Proper exposure using standard surgical procedure.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Good alignment/reduction of the fracture/osteotomy.
- Radiographs can be taken to check the alignment/reduction after fixation.
- Meticulous hemostasis and complete primary soft tissue closure over the implant are essential.

#### Instrumentation

- The **INION FreedomScrew™** products are primarily designed to be used with customized **INION®** instrumentation or general ISO 9714 compatible instrumentation. Corresponding **INION FreedomScrew™** diameters and instrument sizes are presented in the Table 1 (below).

**Table 1. INION FreedomScrew™ diameters and corresponding instrument sizes.**

Screw diameter	K-wire diameter	Drill bit diameter	Tap diameter	Drill tap diameter	Countersink	Screwdriver shaft	Holding sleeve
2.0 mm	-	1.5 mm	2.0 mm	2.0 mm**	for 2.0/2.7 mm screws	hex 2.5 mm	2.0/2.7 mm
2.7 mm	-	2.0 mm	2.7 mm	2.7 mm**			
3.5 mm	-	2.5 mm	3.5 mm	3.5 mm**	for 3.5-4.5 mm screws	hex 3.5 mm	3.5/4.0/4.5 mm
3.5 mm (cannulated)	1.25 mm	2.7 mm (cannulated)	3.5 mm (cannulated)	3.5 mm (cannulated)**			
4.0 mm (cannulated)	1.6 mm	3.2 mm* (cannulated)	4.0 mm (cannulated)	4.0 mm (cannulated)**			
4.5 mm	-	3.2 mm	4.5 mm	4.5 mm**			
4.5 mm (cannulated)	1.6 mm	3.2 mm (cannulated)	4.5 mm (cannulated)	4.5 mm (cannulated)**			

\*INION FreedomScrew™ Ø4.0 mm is designed to be inserted only with 3.2 mm drill (cannulated).

\*\*As customized INION® instruments only

### Fixation with screws

- Choose appropriate INION FreedomScrew™ product for the indication. If needed, the screws can be cut to desired length with the INION® screw cutting pliers. DO NOT use the INION® temperature cautery instrument to cut screws before insertion! DO NOT immerse screws in a hot water bath! **NOTE: Always rinse the drill hole before screw insertion. If needed, use drill bits, bone taps and drill tap instruments through appropriate drill guide. This is to ensure correct hole position and direction, to evaluate appropriate tapping depth and to protect the surrounding tissues from the rotating instrument.**
- Cannulated instruments must be used with cannulated screws. It is recommended to make an x-ray check of the location of K-wire before drilling.
- DRILL AND MANUAL TAP METHOD:
  - Using the appropriate (see Table 1) INION® or other compatible drill bit attached to a slow speed drill (maximum speed of 2000 rpm), drill a suitable hole for the screw into the bone. Use irrigation. If needed, use drill bits through an appropriate drill guide. This is to ensure correct hole position and direction, and to protect the surrounding tissues from the rotating instrument. Take care to hold the drill bit parallel to the long axis of the drill hole while drilling to prevent drill-wobble induced enlargement of the hole.
  - Tap the screw holes manually using the appropriate (see Table 1) INION® or other compatible bone tap before screw insertion. If needed, use bone taps through an appropriate drill guide. This is to ensure correct hole position and direction, to evaluate appropriate tapping depth, and to protect the surrounding tissues from the rotating instrument. Hold the bone tap parallel to the long axis of the drill hole both during insertion and reversing of the bone tap. Only use gentle pressure when tapping to preserve the prepared threads. Care must be taken to tap to a sufficient depth to allow complete seating of the screw. The prepared hole should be irrigated prior to screw insertion to flush out bone debris.
- DRILL TAP METHOD:
  - Instruments for simultaneous drilling and tapping are marked with a yellow stripe on the instrument shaft.
  - Use an appropriate (see Table 1) INION® drill tap instrument attached to a slow speed drill, which has a reverse function, to create a suitable threaded hole for the screw into the bone before screw insertion. Advance the appropriate INION® drill tap slowly using a moderate amount of pressure to push the instrument forward against the bone in the direction parallel to the long axis of the bone hole being created. Use irrigation. Always create the holes for the screws through the drill guide to ensure correct final hole size, and hole position and direction, and to evaluate appropriate tapping depth, and to protect the surrounding tissues from the rotating instrument. Take care to hold the drill tap instrument parallel to the long axis of the bone hole when inserting the screws and when reversing the drill tap instrument during bone hole creation. This prevents hole enlargement induced by instrument wobble, or inappropriate threads on bone. To remove the drill tap from the bone, set the drill machine to reverse first. To preserve the prepared heads, pull gently when reversing the drill tap instrument from the threaded bone hole by using a slow speed drill. Carefully create a sufficient depth for bone hole to allow complete seating of the screw. DO NOT advance drill tap instrument deeper into the bone than its threaded part allows in order to preserve the prepared threads! The prepared hole should be irrigated prior to screw insertion to flush out bone debris.
- SCREW LOADING & PLACEMENT:
  - Attach the appropriate (see Table 1) INION® or other compatible screwdriver shaft with the appropriate holding sleeve. Align the sleeve with an insertion adapter of the screw and connect the sleeve with the insertion adapter by sliding the sleeve towards the screw. Take care to ensure that the sleeve is aligned and connected properly with insertion adapter.
  - Hold the screwdriver shaft and the screw parallel to the long axis of the drill hole and insert the screw fully into the drill hole. DO NOT overtighten! Always insert cannulated screws over a K-wire.
  - When screw is appropriately inserted detach the insertion adapter from the screw head. DO NOT implant the disposable single-use metallic insertion adapter part.
  - When the head of the fully threaded screw is needed as additional support, the appropriate countersink instrument can be used in order to make space for the screw head and to avoid soft tissue irritation by the protruding screw head.
  - When the screw head is not needed (e.g., in the case of syndesmosis fixation or solid cortex fixation), the screw head can be cut by the INION® temperature cautery instrument or surgical cutting pliers to avoid soft tissue irritation. After screw head removal the INION® temperature cautery instrument can be used to finalize and smooth screw shaft remains or to recreate secondary screw head by expanding screw shaft by heat. DO NOT cut or remove the head of a lag screw.
- If a screw gets stuck during its insertion and cannot be fully inserted, remove the screw from the drill hole by reversing the screw out by using the screwdriver. If the screw breaks or the head of the screw gets damaged during screw insertion or reversing and the screw cannot be removed with the screwdriver, drill through the screw with the drill bit that was originally used to create the drill hole for the screw. Tap the drill hole with one size larger (corresponding to the screw diameter) bone tap, and insert a corresponding larger diameter screw into the drill hole following the above described instructions.
- If a screw feels loose, remove the screw from the drill hole by reversing the screw out. Use the screwdriver or drill through the screw with the drill bit that was originally used for creating the drill hole for the screw. Then re-tap the drill hole with one size

larger (corresponding to the screw diameter) bone tap or drill tap, and insert a corresponding larger diameter screw into the drill hole following the above described instructions.

### Fastening of the plates/meshes

- When using **INION FreedomScrew™** products in conjunction with plates or meshes of **INION OTPS™ Systems**, please refer also to the product line-specific instructions for use.
- When **INION FreedomScrew™ Ø2.7/3.5/4.0/4.5 mm** products are used in fastening of the **INION FreedomPlate™**, having drilled and tapped hole through the plate, the head of the screw can be cut off 1-2 mm above the plate surface by the **INION®** temperature cautery instrument or along the plate surface by cutting pliers to provide low-profile plate-screw interface, and therefore to prevent soft tissue irritation. Carefully ensure stability of fixation after screw head is removed. **DO NOT** remove screw head when used with other plate/mesh than **INION FreedomPlate™**. **DO NOT** use **INION FreedomScrew™ Ø2.0 mm** products as a primary fixation screw for the **INION FreedomPlate™**.
- After screw head removal, the **INION®** temperature cautery instrument can be used to finalize and smooth a plate-screw interface and/or to recreate secondary screw head by expanding screw shaft by heat.

### Postoperative reminders

- As with any surgical procedure, careful postoperative management is important for optimal healing.
- Provide the patient with detailed instructions for postoperative care (e.g., regarding immobilization and hygiene maintenance).
- Use appropriate additional immobilization (e.g., a suitable cast, brace and/or crutches) during bone healing.
- Antibiotic therapy at the discretion of the clinician.

### Evaluation of results

Radiographs, magnetic resonance imaging (MRI) and/or computer tomography (CT) scans can be taken postoperatively to evaluate healing. The device is MR Safe.

### WARNINGS

- The **INION FreedomScrew™** products provide fixation and are not intended to replace healthy bone or withstand the stress of full load bearing.
- Incorrect selection, placement, positioning, and fixation of the implant can cause subsequent undesirable results. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.
- Thin soft tissue coverage over the implant may increase the risk of complications.
- Implant may be covered with fibrous tissue capsule when implanted.

### PRECAUTIONS

- Instruments are available to aid accurate implantation of the **INION FreedomScrew™** products. Surgical instruments are subject to wear with normal usage and may break. Surgical instruments are only to be used for their intended purpose. All instruments are to be regularly inspected for wear and damage. Use only the **INION®** or other compatible (see Table 1) instruments.
- **DO NOT** use for unintended applications! Proper function (i.e., effectivity and safety) of these implants cannot be guaranteed in case of off-label use.
- The patient should be warned that the implants can break or loosen as a result of early stress, activity or load bearing. Premature discontinuation of the cast or other immobilization technique may cause non-union or mal-union of the fracture or osteotomy.

### ADVERSE EFFECTS

Complications are similar to those encountered with any other method of internal fixation:

- Premature bending, loosening, fracture or migration of the devices may result from early stress, activity or load bearing.
- Infection can lead to failure of the procedure.
- Neurovascular injuries can occur due to surgical trauma.
- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.

### STERILITY

The **INION FreedomScrew™** products have been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Use only devices that are contained in unopened and undamaged packages. For single-use only. **DO NOT** resterilize or reuse! If the device is re-sterilized or reused, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed. **DO NOT** use implant beyond the expiration date on the label.

### STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F) at a normal relative humidity. Product should not exceed a maximum temperature defined on the indicator. **DO NOT** use if the indicator dot is black.

### CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician.

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
Fax +358 10 830 6601  
info@inion.com  
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

## SYMBOLS USED IN LABELING

	Manufacturer		
	Catalogue number		Lot number
	Date of manufacture		Use by
	Sterilized using irradiation		CE marked product
	Storage temperature range		Keep dry
	Caution, consult accompanying documents		Do not re-use

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INION Oy
Inion® on Inion Oy:n rekisteröity tavaramerkki.