INION CompressOn™ INSTRUCTIONS FOR USE

DESCRIPTION

The **INION CompressOn™** products are cannulated headless compression screws made of degradable co-polymers composed of L-lactic acid and D-lactic acid. These polymers have a long history of safe medical use and they degrade in vivo by hydrolysis into alpha-hydroxy acids that are metabolised by the body. Based on in vitro testing, the screws retain their initial strength up to 12 weeks after implantation and gradually lose their strength thereafter. Bioresorption takes place within two to four years.

The **INION CompressOn™** products are offered in sizes Ø 2.7–6.0 mm and are designed to be used with **INION®** instrumentation: bone drill bits, drill taps, screwdriver shafts, drill sleeves, K-wires and depth gauges. All screws can be inserted over a guide wire.

The **INION CompressOn™** products are sterile, non-pyrogenic and non-collagenous.

Product selection

INION CompressOn™ Screw and corresponding **INION®** instrument sizes

Screw size	Nonsterile, Reusable Instruments				
	K-wire size	Drill sleeve	Cannulated Drill bit	Cannulated Drill tap	Screwdriver shaft
2.7 mm	0.8 mm	IFS-9022	IFS-2703 (Ø 2.0 mm)	IFS-2708 (Ø 2.7 mm)	IFS-9003
3.5 mm	1.25 mm	IFS-9023	IFS-3503 (Ø 2.7 mm)	IFS-3508 (Ø 3.5 mm)	IFS-9004
4.0 mm	1.25 mm	IFS-9024	IFS-4513 IFS-4503 (Ø 3.2 mm)	IFS-4018 IFS-4008 (Ø 4.0 mm)	IFS-9004
4.5 mm	1.6 mm	IFS-9024	IFS-4513 IFS-4503 (Ø 3.2 mm)	IFS-4518 IFS-4508 (Ø 4.5 mm)	IFS-9005
5.0 mm	1.6 mm	IFS-9025	IFS-5003 (Ø 3.7 mm)	IFS-5008 (Ø 5.0 mm)	IFS-9005
6.0 mm	2.0 mm	IFS-9026	IFS-6003 (Ø 4.2 mm)	IFS-6008 (Ø 6.0 mm)	IFS-9006

INDICATIONS

The **INION CompressOn™** products are intended for maintenance of alignment and fixation of bone fractures, comminuted fractures, osteotomies, arthrodeses or bone grafts (i.e., autografts or allografts) in the presence of appropriate additional immobilization (e.g., rigid fixation implants, cast or brace).

In addition, the **INION CompressOn™** Ø 3.5/4.0/4.5/5.0/6.0 mm products are specifically intended for use in the following indications:

- A. General indications: maintenance of reduction and fixation of cancellous bone fractures, osteotomies or arthrodeses of the upper extremity, ankle and foot in the presence of appropriate brace and/or immobilization.
- B. Specific indications: fractures and osteotomies of the malleoli, and ankle fractures.

CONTRAINDICATIONS

The **INION CompressOn™** products are not intended for use in and are contraindicated for:

- Active or potential infection.
- Patient conditions including limited blood supply, insufficient quantity or quality of bone; and where patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse).
- Spine indications.
- High-load bearing indications (e.g., diaphyseal fractures of long bones) unless used in conjunction with traditional rigid fixation.

INFORMATION FOR USE

Surgical considerations and reminders

- Prophylactic perioperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.
- Maintain sterile field throughout the procedure.
- Proper exposure using standard surgical procedure.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Good alignment/reduction of the fracture/osteotomy.

- Radiographs can be taken to check the alignment/reduction after fixation.
- Meticulous hemostasis and complete primary soft tissue closure over the implant are essential.

Fixation with screws

PREPARATION OF A THREADED BONE HOLE:

1. Select appropriate **INION CompressOn™** product/diameter for the indication.
 2. Use one or multiple K-wires to reduce the bone fragments. It is recommended to make an X-ray check of the location of K-wire before proceeding. Measure the needed bone hole depth using appropriate depth gauge over a K-wire and select correct screw length accordingly. Ensure that selected screw can be inserted completely in to the final compressed bone hole.
 3. Use drill tap instrument through appropriate drill sleeve. This is to ensure correct hole position and direction, to evaluate appropriate tapping depth and to protect the surrounding tissues from the rotating instrument.
 4. Use an appropriate cannulated **INION®** drill tap instrument attached to a slow-speed power drill, which has a reverse function to create appropriate threaded bone hole. Instruments for simultaneous drilling and tapping (drill taps) are marked with a yellow stripe on the instrument shaft.
 5. Advance the drill tap over a K-wire slowly using moderate pressure. Push the instrument forward by maintaining the axial alignment. Use irrigation. **DO NOT** advance drill tap instrument deeper than its threaded part.
 6. Evaluate appropriate bone hole depth by reading the depth markings on drill tap. Ensure the bone hole is long enough to accommodate the selected screw in the final compressed bone hole.
 7. To remove the drill tap from the bone, set the slow speed power drill to reverse. Maintain the axial alignment when reversing the drill tap instrument to prevent enlarged bone hole or inappropriate threads.
- **NOTE:** For drilling and tapping separately, follow steps 1 and 2. In step 3 use **INION®** drill bit to create pilot hole. Tap the pilot hole using compatible cannulated **INION®** drill tap before screw insertion. Maintain axial alignment in all steps to prevent enlarged bone hole or inappropriate threads.

SCREW LOADING & PLACEMENT:

8. Connect appropriate screwdriver shaft with the handle.
 9. Attach the selected screw to the corresponding **INION®** screwdriver shaft.
 10. Always rinse the drill hole before screw insertion to remove any debris.
 11. Insert the screw fully into the bone hole by maintaining axial alignment in order to avoid soft tissue irritation by protruding screw head. Apply moderate pressure to the screw with screwdriver while screw head is advancing in to the bone hole. **DO NOT** use excessive force. **DO NOT** advance the screw head past the fracture/osteotomy line. Always insert screws over a K-wire.
 12. In case of partly protruding screw head after insertion, shape the protruding part as needed to avoid soft tissue irritation, by using e.g. pliers, oscillating saw or a **INION®** temperature cautery instrument.
 13. When screw is appropriately inserted detach the screwdriver from the screw by pulling it out axially to avoid damaging the screw head or screwdriver.
- **NOTE:** If screw gets stuck during insertion and cannot be fully inserted, remove the screw by reversing using the screwdriver. If screw gets damaged during insertion and cannot be removed by reversing and fixation feels unstable, drill through the screw with the drill bit which corresponds to damaged screw or one size larger. Tap the drill hole with one size larger (corresponding to the screw size) drill tap, and insert a corresponding one size larger diameter screw into the drill hole following the above described instructions.

Postoperative reminders

- As with any surgical procedure, careful postoperative management is important for optimal healing.
- Provide the patient with detailed instructions for postoperative care (e.g., regarding immobilization and hygiene maintenance).
- Use appropriate additional immobilization (e.g., a suitable cast, brace and/or crutches) during bone healing.
- Antibiotic therapy at the discretion of the clinician.

Evaluation of results

Radiographs, magnetic resonance imaging (MRI) and/or computer tomography (CT) scans can be taken postoperatively to evaluate healing. **INION CompressOn™** Screw is MR Safe.

WARNINGS

- The **INION CompressOn™** products provide fixation and are not intended to replace healthy bone or withstand the stress of full load bearing.
- Incorrect selection, placement, positioning, and fixation of the implant can cause subsequent undesirable results. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.
- Do not cut the screw before insertion. When necessary, e.g. pliers, oscillating saw or a temperature cautery instrument can be used to shape partly protruding screw head of **INION CompressOn™** screw after insertion. In such a case, shape the protruding part as needed to avoid soft tissue irritation.
- Do not resterilize the implant.
- Do not reuse the implant.
- Do not use implant beyond expiration date printed on the label.
- Do not use devices that are contained in opened or damaged packages.
- **USE ONLY IF** sterilization indicator shows exposed color (red).
- Remove the implant from the package using aseptic methods to retain implant sterility.

PRECAUTIONS

- Instruments are available to aid accurate implantation of the **INION CompressOn™** products. Surgical instruments are subject to wear with normal usage and may break. Surgical instruments are only to be used for their intended purpose. All instruments are to be regularly inspected for wear and damage. Use only **INION®** instruments (see Product Selection table).
- **DO NOT** use for unintended applications! Proper function (i.e., effectivity and safety) of these implants cannot be guaranteed in case of off-label use.
- The patient should be warned that the implants can break or loosen as a result of early stress, activity or load bearing. Premature discontinuation of the cast or other immobilization technique may cause non-union or mal-union of the fracture or osteotomy.

ADVERSE EFFECTS

Complications are similar to those encountered with any other method of internal fixation:

- Premature bending, loosening, fracture or migration of the devices may result from early stress, activity or load bearing.
- Infection can lead to failure of the procedure.
- Neurovascular injuries can occur due to surgical trauma.
- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.
- Implant may be covered with fibrous tissue capsule when implanted.

STERILITY

The **INION CompressOn™** products have been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Use only devices that are contained in unopened and undamaged packages. For single-use only. **DO NOT** re-sterilize or re-use! If the device is re-sterilized or reused, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed.

STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F) at a normal relative humidity. Product should not exceed a maximum temperature defined on the indicator. **DO NOT** use if the indicator dot is black.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician.

MANUFACTURER

INION Oy
Lääkärintäti 2
33520 Tampere
Finland
Tel. +358 10 830 6600
Fax +358 10 830 6601
info@inion.com
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

	Manufacturer		Prescription Use Only
	Catalogue number		Lot number
	Date of manufacture		Use by
	Sterilized using irradiation		CE marked product
	Storage temperature range		Keep dry
	Caution, consult accompanying documents		Do not re-use
	Do not re-sterilize		Consult instructions for use

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Inion®, Inion Oy'un tescilli ticari markasıdır.