ENGLISH

INSTRUCTIONS FOR USE INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM

DESCRIPTION

The **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM** consists of plates and screws made of degradable co-polymers composed of L-lactic and D-lactic acid. These polymers have a long history of safe medical use and they degrade in vivo by hydrolysis into alpha-hydroxy acids that are metabolised by the body. Based on in vitro testing, the implants retain most of their initial strength up to 16 weeks and gradually lose their strength thereafter; and bioresorption takes place within two to four years. The plates SPN-5110 – SPN-5121 and screws SPN-5320 – SPN-5325 contain radiopaque tantalum markers for postoperative radiographic imaging.

The **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM** implants are available in different sizes and are designed for use with customised instrumentation, e.g., the **INION S-1™** drill bits, threading taps, awl, plate holder, drill guides, screw driver, temporary fixation screw, screw removal instrument, plate templates, and a heating device.

The **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM** implants are sterile and non-collagenous.

INDICATION

The **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM**, in conjunction with adjunctive immobilisation (i.e., simultaneous conservative treatment with a collar), is intended for use in anterior cervical discectomy and fusion procedures as a means to maintain the relative position of block bone allografts or autografts. This device is not intended for use in other areas of the spine, and is not for load bearing indications.

CONTRAINDICATIONS

The **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM** implants are not intended for use in, and are contraindicated for:

- Active or potential infection
- Cancer
- Pseudoarthrosis
- Patient conditions including limited blood supply, insufficient quantity or quality of bone; and where patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse)

INFORMATION FOR USE

Surgical considerations and reminders

- Prophylactic preoperative antibiotic treatment is recommended.
- Maintain a sterile field throughout the procedure.
- Proper exposure using standard surgical procedure for anterior cervical discectomy and fusion.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Good approximation of tissues.
- Meticulous hemostasis and complete primary soft tissue closure over the implant are essential.

Preparation of the plate

A. Before starting the procedure (see also the instructions for use of the **INION THERMO+™** water bath):

- Place a sterile **INION THERMO™** drape to the **INION THERMO+™** water bath. DO NOT switch on the water bath until filled with liquid!
- Pour sterile water or saline (min 0.4 litres max 0.7 litres) into the waterbath.
- Switch the water bath on. Wait until the liquid heats up to the setting temperature 70°C / 158°F; this takes about 30 minutes.

B. Following anterior cervical discectomy and fusion procedure with appropriate block bone graft:

- Use the **INION S-1™** plate templates to determine the appropriate plate size. Place the appropriate sized template against bone at the intended location of the plate to confirm correct plate length and screw hole position relative to the endplates. The plate should be positioned to cover the bone graft to maintain the relative position of the graft. DO NOT implant the template!
- After the appropriate **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM** plate has been determined, the plate should be contoured to the bone. Immerse the plate in the **INION THERMO+™** water bath for approximately one minute. After water bath treatment, the plate can be contoured to the bone also by manually pressing it directly against the bone at its intended location. The transparency of the plate helps in the determination of appropriate positioning. Contour the plate to desired shape by hand and/or by using pliers. The plate is malleable for approximately 5-10 seconds and can be shaped as desired. The plate must be contoured to restore and maintain the three-dimensional anatomical alignment and balance of the spine. Please note that the plate must always be implanted with the countersunk side upwards!
- Maximum total immersion time of the plate/mesh in the waterbath is 30 minutes. Prolonged immersion may weaken the plate/mesh.
- Ensure that the plate is accurately adapted to the bone contours.
- If contouring is not satisfactory, the water bath treatment can be repeated by immersing the plate in the water bath for additional 20-30 seconds. If only part of the plate needs further adaptation, dip only that portion of the plate.
- If the wound is too deep to allow the plate to be contoured directly against the bone by hand or if the plate is not accurately adapted to the spine after contouring, the plate template can be used to contour the plate. The template should first be contoured to the bone. Contour the template with pliers and confirm the correct final shape by placing the contoured template against bone at the intended

location of the plate. Ensure that the contoured template is finally adapted accurately to the bone contours, and remove the template thereafter. Place the plate over the contoured template and keep them gently pressed together with pliers. Immerse the plate and the contoured template together in the **INION THERMO+**™ water bath keeping the plate and the template over each other with pliers. Take care not to apply too much force on the plate with the pliers; this is to avoid creating depressions on the plate. After water bath treatment the plate can be formed to the shape of the contoured template by pressing the plate against the template by hand and by using the **INION**® plate bending pliers. Place the contoured plate against the bone at its intended location and ensure that the plate is accurately adapted to the bone contours.

- The **INION S-1**™ plate holder instrument can be used to hold both the plates and plate templates. The plate holder instrument can be connected to the plate or template either before or after the contouring of the plate or template. This is done by inserting the tip of the instrument into the middle plate-hole or template-hole intended for the plate holder and tightening the screw on the plate holder.
- DO NOT bend the plate before or without water bath treatment, or after the plate has cooled! Repeated bending can weaken the plate, especially once it has cooled. Bending while warm can be carried out up to three times.
- Take care not to perforate the **INION THERMO**™ drape with surgical instruments.
- DO NOT use implants that show signs of damage!
- DO NOT immerse the screws in the water bath!

Fastening of the plate

- After contouring, the plate should be held in place with the plate holder. Additionally, insert at least two temporary fixation screws through opposite screw holes of the plate to hold it in place.
- Insertion of the temporary fixation screw:
 1. Make sure the plate is in the correct position. Ensure that the countersunk side of the plate is upwards!
 2. Place the tip of the **INION S-1**™ drill guide into the intended screw hole of the plate, angulate the drill guide in the desired direction for the screw, and maintain the position and direction of the drill guide thereafter.
 3. Before inserting the temporary screw, make a starting hole through the cortex of the bone by using the **INION S-1**™ awl instrument. Pass the awl instrument through the drill guide and the hole of the plate. This is to ensure correct hole position and direction, and to protect the surrounding tissues from the awl. Maintain the position and direction of the drill guide and gently tap the awl with a small mallet to create a starting hole. NOTE: Always tap the awl only as long as the marking line on the awl's shaft remains above the top of the drill guide. Tapping any further may result in plate breakage by the awl.
 4. Align the screwdriver of the metallic temporary screw with the top of the screw. Push the screw head into the corresponding socket of the screwdriver. Insert the self-tapping temporary fixation screw through the drill guide and the hole of the plate, placing it fully into the drill hole. While doing this it is important to still maintain the position and direction of the drill guide; hold the screwdriver and the screw parallel to the long axis of the drill hole and keep the plate in position with the plate holder. DO NOT overtighten!
- After inserting the temporary fixation screws, loosen the tip of the plate holder instrument from the hole of the plate and remove the plate holder.
- Choose appropriate **INION S-1**™ 4.5-mm screws for monocortical fixation of the plate. Preoperative X-rays can be used to determine appropriate screw length for monocortical fixation but the magnifying effect of X-rays should be taken into account when estimating the screw length (i.e., choose screws slightly shorter than indicated by the X-ray based measurements).
- Place the tip of the **INION S-1**™ drill guide into the intended screw hole of the plate, angulate the drill guide in the direction desired for the screw, and maintain the position and direction of the drill guide thereafter. However, take care not to angulate the drill guide too much to avoid drilling through the other screws. Please note that if a previously inserted screw is hit with the drill when drilling holes for additional screws, the radiopaque marker located in the tip of the screw may become displaced or removed. Angulation of the screw may also cause protrusion of the screw head.
- Before drilling the hole:
 1. Make a starting hole for the drill bit through the cortex of the bone by using the **INION S-1**™ awl instrument. Gently tap the awl instrument through the drill guide and the hole of the plate with a small mallet - to create a starting hole - and thereafter remove the awl instrument but maintain the position and direction of the drill guide. Always apply the awl instrument through the drill guide and the hole of the plate to ensure correct hole position and direction, and to protect the surrounding tissues from the awl.
 2. Choose the appropriate screw-specific **INION S-1**™ drill bit and bone tap instruments with stops (corresponding to the screw diameter and length). Please note the colour and stripe codes on the instrument shafts:

Green	Screw length 12 mm
Blue	Screw length 14 mm
Magenta	Screw length 16 mm
One colour stripe	Screw diameter 4.5 mm
Two colour stripes	Screw diameter 5.0 mm
- Using the appropriate screw-specific **INION S-1**™ drill bit with a stop (corresponding to the screw diameter and length) attached to a slow speed drill (maximum speed of 2000 rpm), complete the monocortical hole for the screw by drilling through the drill guide and the hole of the plate until the stop of the drill bit reaches the drill guide. Use irrigation. Always drill the holes for the screws through the drill guide and the holes of the plate to ensure correct final plate-hole size, drill-hole position, direction and depth, and to protect the surrounding tissues from the drill. Always use the screw-specific **INION S-1**™ drill bit with the stop (corresponding to the screw diameter and length) to avoid drilling too deep. Care must be taken to drill until the stop reaches the drill guide; this allows complete seating of the screw. Take care to hold the drill guide and the bone drill parallel to the long axis of the drill hole while drilling to prevent enlargement of the hole induced by drill-wobble.
- Using the appropriate screw-specific **INION S-1**™ bone tap with a stop (corresponding to the screw diameter and length), insert the bone tap into the screw hole manually through the drill guide and the hole of the plate until the stop of the bone tap reaches the drill guide. Always tap the screw hole before screw insertion, and always insert the bone tap into the screw hole through the drill guide and the hole of the plate. This ensures the tapping of the plate-hole, the correct tapping position, direction and depth, and protects the surrounding tissues from the bone tap. During bone tap insertion, press the plate tightly against bone either by hand or by using the drill guide and the **INION S-1**™ plate holder instrument. Always use the screw-specific **INION S-1**™ bone tap with the stop (corresponding to the screw diameter and length) to avoid tapping too deep. Care must be taken to insert the bone tap until the stop reaches the drill guide to allow complete seating of the screw. However, to preserve the prepared threads DO NOT continue insertion of the bone tap when the stop has reached the drill guide. The position and depth of the drill hole can be confirmed by imaging before reversing the bone tap when the bone tap is fully advanced into the drill hole (i.e., when the stop of the bone tap reaches the drill guide). Take care to hold the

drill guide and the bone tap parallel to the long axis of the drill hole both during insertion and reversing of the bone tap, and use only gentle pressure when tapping to preserve the prepared threads.

- Rinse the drill hole before screw insertion to flush out bone debris and to enable complication-free screw insertion.
- Align the **INION S-1™** screwdriver with the top of the previously selected screw. Gently push the screwdriver tip into the corresponding socket of the screw head. Ensure that the screwdriver tip is fully inserted into the screw socket by pressing it firmly against the screw. Ensure that the screwdriver is aligned properly before pushing it into the screw, since there is a risk of damaging the plastic screw head. Still maintaining the position and direction of the drill guide, hold the screwdriver and the screw parallel to the long axis of the drill hole and while keeping the plate in position and tightly against bone surface, insert the screw through the drill guide and the hole of the plate fully into the drill hole. DO NOT overtighten!
- Create additional threaded drill holes and insert additional screws through the remaining screw holes of the plate following the above-described instructions for screw insertion. Always insert screws through all screw holes of the plate! However, DO NOT insert screws through the middle plate-hole(s) intended for the plate holder instrument! The temporary fixation screws should be removed before creating the drill hole for the last screws. Never leave the temporary fixation screws in situ postoperatively.
- If a screw feels loose, remove the screw, tap the drill hole by using the appropriate **INION S-1™** 5.0-mm bone tap (corresponding to the length of the original screw), and insert the appropriate **INION S-1™** 5.0 mm emergency screw (corresponding to the length of the original screw) into the drill hole following the above-described instructions for screw insertion. Please note that the 5.0 mm screw must be inserted without the drill guide.
- If a screw gets stuck during its insertion and cannot be fully inserted, remove the screw from the drill hole by reversing the screw out using the screw driver. If the screw breaks or the head of the screw gets damaged during screw insertion or reversing and the screw cannot be removed with the screw driver, use the **INION S-1™** screw removal instrument (see instructions below), or drill through the screw with the same screw-specific drill bit that was originally used to create the drill hole for the screw. Always use the drill guide when drilling through screws and take care to angulate the drill guide and the drill parallel to the long axis of the screw. Please note that when a screw is drilled out, the tantalum marker located in the tip of the screw may become displaced and/or remain in the bone. This should be kept in mind when reading postoperative radiographs.
- Screw removal with the screw removal instrument:
 1. Using the **INION S-1™** 2.0 mm drill bit attached to a slow speed drill (maximum speed of 2000 rpm), drill a hole for the screw removal instrument centrally into the screw intended to be removed. Drill through the drill guide and the hole of the plate until the stop of the drill bit reaches the surface of the plate or head of the screw. Use irrigation. Always drill the hole for the screw removal instrument through the drill guide and the hole of the plate to ensure correct drill-hole position and direction, and to protect the surrounding tissues from the drill. Always use the **INION S-1™** 2.0 mm drill bit with the stop to avoid drilling too deep.
 2. Withdraw the outer sleeve of the **INION S-1™** screw removal instrument by turning it clockwise until it touches the instrument handle. Insert the tapered end of the **INION S-1™** screw removal instrument into the drill hole by turning counter-clockwise. Keep turning until the instrument reaches the bottom of the drill hole. Then turn the outer sleeve counter-clockwise until it reaches the surface of the screw. Thereafter turn the instrument handle with the sleeve counter-clockwise to remove the biodegradable screw. Always use the screw removal instrument with the drill guide to protect the surrounding tissues from the screw removal instrument.
- After screw removal either by reversing the screw with the screw driver or the screw removal instrument:
 1. Rinse the drill hole, and re-tap the drill hole with the correct screw-specific **INION S-1™** bone tap (corresponding to the screw diameter and length) following the above-described instructions. Care must be taken to insert the bone tap until the stop of the bone tap reaches the drill guide to allow complete seating of the screw. However, DO NOT continue insertion of the bone tap after the stop has reached the drill guide to preserve the prepared threads. If the bone tap cannot be inserted deep enough (i.e., until the stop reaches the drill guide), choose a shorter screw or remove the tap and re-drill the hole with the correct screw-specific drill bit (corresponding to the screw diameter and length) following the above-described instructions. Care must be taken to drill until the stop of the drill reaches the drill guide; this allows complete seating of the screw. Re-tap the hole with the correct screw-specific bone tap. Always use the drill guide when re-drilling and re-tapping, and take care to angulate the drill guide and the instruments parallel to the long axis of the original drill hole. Re-drilling or re-tapping in a different angle than that chosen for the original hole may damage the previously created threads and cause enlargement of the hole.
 2. Rinse the drill hole before screw insertion to flush out bone debris and to enable complication-free screw insertion.
 3. Insert the appropriate **INION S-1™** screw following the above-described instructions for screw insertion. The previously removed screw can be used if it is undamaged.
- After screw removal by drilling through the biodegradable screw:
 1. Rinse the drill hole and tap it with the appropriate screw-specific **INION S-1™** 5.0-mm bone tap (corresponding to the original screw length) following the above-described instructions. Always use the drill guide when drilling and tapping, and take care to angulate the drill guide and the instruments parallel to the long axis of the original screw and drill hole. Drilling or tapping in a different angle than that chosen for the original hole may cause enlargement of the hole. If the drill hole feels too tight for the 5.0-mm bone tap, re-drill the hole first with the appropriate **INION S-1™** 4.2-mm drill bit (corresponding to the original screw length).
 2. Insert the appropriate **INION S-1™** 5.0-mm emergency screw (corresponding to the original screw length) following the above-described instructions for screw insertion.

Postoperative reminders

- As with any surgical procedure, careful postoperative management is important for optimal healing.
- Provide the patient with detailed instructions for postoperative care (e.g., regarding appropriate immobilisation and hygiene maintenance).
- Use appropriate additional immobilisation (i.e., a suitable collar) until bone healing has been observed in radiographic examination. Typical immobilization time is 3-6 weeks.
- Antibiotic therapy at the discretion of the clinician.

Evaluation of results

- Radiographs, magnetic resonance imaging (MRI) and/or computer tomography (CT) scans can be taken postoperatively to evaluate healing.
- The plates SPN-5110 – SPN-5121 and screws SPN-5320 – SPN-5325 contain radiopaque tantalum markers to enable radiographic evaluation of implant position. The plates contain one marker on both ends of the plate (i.e., totally two markers per plate) and screws have a marker in the tip of the screw. Lateral radiographs are most suitable for evaluation of marker/implant position. Please note that if

a screw has been removed by drilling it out, the marker of the removed screw may have got displaced and/or remain in the bone. Please also note that screw markers may become displaced or removed if a previously inserted screw is hit with the drill when drilling holes for additional screws. When evaluating position of tips of the screws, please note that the created drill hole is always slightly longer than the corresponding screw. When evaluating implant position, always compare the radiographic images to the x-rays taken immediately postoperatively. If implant displacement or breakage is suspected based on the location of radiopaque markers seen on radiographs, always confirm the finding with MRI and/or CT before considering reoperation!

WARNINGS

- The **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM** devices provide fixation and are not intended to replace healthy tissues or withstand the stress of full load bearing.
- Incorrect selection, placement, positioning, or fixation of the implant can cause subsequent undesirable results or breakage of implants or instruments. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.

PRECAUTIONS

- Instruments are available to aid accurate implantation of the **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM**. Surgical instruments are subject to wear with normal usage and may break. Surgical instruments are only to be used for their intended purpose. All instruments are to be regularly inspected for wear and damage. Use only the **INION S-1™** instruments.
- Use only **INION S-1™** screws to secure the **INION S-1™** plates. DO NOT use **INION S-1™** screws to secure plates or sheets other than those of **INION®** (e.g., metal plates or sheets).
- DO NOT use for unintended applications! The proper function (i.e., efficacy and safety) of these implants cannot be guaranteed in the event of off-label use.
- The patient should be warned that the implants can break or loosen as a result of early stress, activity or load bearing. Premature discontinuation of additional external immobilisation (i.e., collar) may cause non-union or mal-union.
- The surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.

ADVERSE EFFECTS

Complications are similar to those with any other method of internal fixation:

- Premature bending, loosening, breakage or migration of the devices may result from early stress, activity or load bearing.
- Infection can lead to failure of procedure.
- Neurovascular injuries can occur due to surgical trauma.
- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances.

STERILITY

The **INION S-1™ BIODEGRADABLE ANTERIOR CERVICAL FUSION SYSTEM** implants have been sterilised with ionising irradiation. Use immediately after opening the sterile seal. Use only devices that are contained in unopened and undamaged packages. For single use only. **DO NOT re-sterilize or reuse**. If the device is re-sterilized or reused, the properties of the implant will be affected, and the safety and efficacy cannot be guaranteed. DO NOT use implant beyond the expiration date on the label.

STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F) at a normal relative humidity. Product should not exceed a maximum temperature defined on the indicator. Do not use if the indicator dot is black.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician.

MANUFACTURER

INION Oy
Lääkärintäti 2
33520 Tampere
Finland
Tel. +358 10 830 6600
Fax +358 10 830 6601
info@inion.com
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

Manufacturer

Catalogue number

Date of manufacture

Sterilized using irradiation

Storage temperature range

Caution, consult accompanying documents

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Inion® é uma marca comercial registada da Inion Oy.