ENGLISH

INSTRUCTIONS FOR CLEANING, DISINFECTION AND STERILIZATION OF REUSABLE INION® INSTRUMENTATION

Cleaning, disinfection and sterilization are required before the first use of non-sterile **INION® INSTRUMENTATION** and before each re-use of all reusable **INION® INSTRUMENTATION**. Effective cleaning and disinfection are a vital requirement for an effective sterilization of the **INION® INSTRUMENTATION**. Cleaning, disinfection and sterilization of reusable **INION® INSTRUMENTATION** must be conducted according to the procedures outlined in this document.

HOW SUPPLIED

INION® INSTRUMENTATION is supplied sterile or non-sterile. Non-sterile **INION® INSTRUMENTATION** must be cleaned and sterilized prior to use according to the procedures outlined in this document.

WARNINGS

- Do not expose the **INION® INSTRUMENTATION** to temperatures higher than 140 °C (284 °F).
- Do not use high **acidic or high alkaline solutions**.
- Do not use saline solutions.
- Do not use any instrument as a lever arm.
- Use the instrument only for its intended use.
- Do not change limb position while instrument is in the joint.
- Do not penetrate soft tissue with the suture grasper.
- Suture passer, 45°: Moving the instrument back and forth through the tissue may result in tip breakage or enlarged tear in the soft tissue. Use for tissue penetration only with instrument jaws closed.
- Excess load directed to the instruments may cause breakage of the instrument tip.
- Careless use of instruments with sharp edges may result in soft tissue damage.

LIMITATIONS ON REPROCESSING

Repeated processing according to the instructions has a minimal effect on the **INION® INSTRUMENTATION**. Life cycle is normally determined by wear and damage due to use. Damaged **INION® INSTRUMENTATION** must not be further used and should be disposed of.

INSTRUCTIONS: POINT OF USE

- Do not allow blood or soil to dry on the instruments.
- Remove excess soil with disposable, non-shedding wipe.
- Flush cannulated devices with cool water.
- It is recommended that instruments are reprocessed as soon as is reasonably practical following use.

MANUAL CLEANING/DISINFECTION

CLEANING:

Equipment: Soaking bath, detergent, brush, running water/purified water, ultrasonic washer unit, DI water

Ensure the compatibility of the cleaning agent with the **INION® INSTRUMENTATION**.

1. Instruments with removable parts must be disassembled prior to cleaning.
2. Immerse the instrument in a non-foaming, neutral enzymatic cleaner or neutral detergent (7.0-9.0 pH) solution and allow to soak for a minimum of 10 minutes. Follow the cleaner or detergent manufacturer's instructions for use for correct exposure time, temperature, water quality and concentration. Alkaline agents with a pH of 12 or less may be used to clean stainless steel and polymer instruments in countries where required by law or local ordinance; or where prion diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeld-Jakob Disease (CJD) are a concern. The capacity of a detergent to inactivate prions is most likely with pH-values > 10 and a soaking time of more than 10 minutes under elevated, but not protein-fixing temperatures (for example 55 °C). However, claims that a product can inactivate prions must be evidence-based ¹⁾. This treatment does not have adverse effect to stainless steel instruments. Aluminium and plastic parts may suffer visual deterioration.
3. Clean the instrument with a soft-bristled brush and common hospital cleaning tools for a minimum of five minutes. Pay extra attention to challenging design features and areas, where tissue residue accumulates the most easily. Move and/or retract all moveable device features and use brush to clean hinges and mating surfaces. Scrub lumens and holes with a brush. Clean long, narrow cannulations meticulously with a brush or by injecting water through the cannula.
4. Rinse under clean running water for min. 2 minutes. Ensure that running water passes through cannulations and that blind holes are repeatedly filled and emptied. Movable features should be actuated during rinsing. It is critical that alkaline cleaning agents are completely and thoroughly neutralized and rinsed from the devices or degradation may occur that limits the device life.
5. Sonicate the instrument in warm enzymatic detergent in its fully opened position for a minimum of 10 minutes.
6. Rinse the instrument thoroughly with warm deionized/purified water for a minimum of 3 minutes. Ensure that running water passes through cannulations and that blind holes are repeatedly filled and emptied. Movable features should be actuated during rinsing.

¹⁾Hygiene Requirements for the Reprocessing of Medical Devices, Robert Koch Institute and BfArM. Bundesgesundheitsbl 2012 • 55:1244–1310.


DISINFECTION:

Equipment: Soaking bath, disinfectant, syringes min. 50 ml, deionized water

Ensure the compatibility of the disinfectant with the **INION® INSTRUMENTATION**.

1. Soak the disassembled products for a minimum of 10 minutes in the disinfectant solution so that the products are sufficiently covered. Pay attention that there is no contact between the products. Sway movable parts several times during disinfection.
2. Rinse all lumens of the products at least five times at the beginning of the soaking time by application of a single-use syringe (minimum volume 50 ml).
3. Remove the products of the disinfectant solution and post-rinse them at least five times intensively (at least 1 min) with deionized water.
4. Rinse all lumens of the products at least five times by application of a single-use syringe (minimum volume 50 ml).
5. Dry the instruments thoroughly by clean, soft, lint-free cloth. Move and/or retract all movable parts during drying. Clean compressed air may be used for drying of inaccessible features. If additional drying is required, place the instruments in a clean area or heat in an oven below 90°C (194 °F).

Figure 1: Injecting water through the cannula (Inion Shoulder instruments)

AUTOMATED CLEANING/DISINFECTION

Automatic washer/disinfector using thermal disinfection is recommended for **INION® INSTRUMENTATION**. Follow the internal hospital instructions and recommendations of manufacturers of automatic washer/disinfector, especially in regard to operating time and temperature. Equipment: Automatic washer/disinfector, cannula brush, running water, cleaner, detergent, deionized or purified water.

Ensure the compatibility of the cleaning agent with the **INION® INSTRUMENTATION**.

- Do not expose to temperatures higher than 140 °C (284 °F).
- INION® INSTRUMENTATION** is classified as critical according to the Ao classification as defined in ISO 15883 Washer-disinfectors.
- Sterile / deionized / purified water is recommended for final rinsing.
 - Conduct steps 1-6 from section Manual Cleaning above.
 - Place the instrument in the automatic washer/disinfector so that the instrument cannulations, hinges, lumens and holes can drain and use a washing cycle with the following minimum parameters:

Pre-wash	5 minutes	Cold drinking water, <43 °C (109 °F)
Wash I	5 minutes	Neutral pH (7.0-9.0) enzymatic cleaner, warm drinking water (>45 °C / 113 °F)
Wash II	5 minutes	Neutral pH (7.0-9.0) detergent wash, warm drinking water (>45 °C / 113 °F)
Rinse	5 minutes	Warm DI or PURW (>40 °C / 104 °F)
Thermal disinfection	5 minutes	90°C (194 °F) ²⁾
Drying	40 minutes	

²⁾ Providing the medical instrument is heat resistant, an autoclave cycle at 134 °C for 18 minutes or more is generally effective for prion inactivation (Fichet et al. 2004: Novel methods for disinfection of prion-contaminated medical devices. Lancet 364, 521–6.)

INSPECTION

- Instruments should be carefully inspected prior to sterilization in order to ensure that all visible contamination has been removed.

REASSEMBLY

- If appropriate, reassemble instruments before sterilization.

FUNCTIONAL TESTING

- Actuate the moving parts to ensure smooth operation throughout the intended range of motion.
- If the instrument is part of a larger assembly, check that devices can be assembled with mating components.
- Examine the instruments for signs of wear, corrosion and non-functioning parts to ensure proper performance of the device.

STERILE PACKAGING

- Packaging should ensure sterility of the instruments until opened for use at the sterile field.
- The **INION® INSTRUMENTATION** must be packaged separately into a single-use sterilization packaging suitable for steam sterilization or in the **INION®** sterilization tray designed for the purpose.

STERILIZATION

- Steam sterilization with pre-vacuum procedure is recommended for the **INION® INSTRUMENTATION**.
- The **INION® INSTRUMENTATION** has to be cleaned and disinfected prior to sterilization.

The recommended steam sterilization parameters are as follows:

Sterilization method	Temperature	Exposure time	Drying time in wrap tray cover off
Pre-vacuum	121 °C 250 °F	20 minutes	60 min
Pre-vacuum	132 °C 270 °F	3 minutes	60 min
Pre-vacuum	134 °C 273 °F	3 minutes	60 min

STORAGE

- Store sterilized instruments in a cool, dry place protected from dust and direct sunlight.
- Examine carefully prior to opening the package that the package integrity has not been compromised.

These instructions have been validated by Inion Oy as being capable of preparing the **INION® INSTRUMENTATION** for re-use. The end user is responsible for the sterility of the **INION® INSTRUMENTATION**. Any deviation from these instructions should be properly evaluated for effectiveness and potential adverse consequences. Follow the internal hospital instructions and recommendations by manufacturers of detergents and of automatic washers and sterilizers.

MANUFACTURER

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33520 Tampere
Finland
Tel. +358 10 830 6600
Fax +358 10 830 6601
info@inion.com
www.inion.com

CE / CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

Medical device

Unique Device Identifier

Manufacturer

Non-sterile product

Consult instructions for use

Lot number

Catalogue number

Date of manufacture

Keep dry

Caution, consult accompanying documents

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INION Oy
Inion® Inion Oy'un tescilli markasıdır.