ENGLISH

INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM INSTRUCTIONS FOR USE

DESCRIPTION

The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** is a crano-maxillofacial plating system intended to maintain accurate alignment of bone fractures and osteotomies in the cranium, mid-face and mandible in conjunction with appropriate maxillomandibular fixation.

The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** consists of fixation plates, meshes and screws made of degradable co-polymers composed of L-lactic acid, D-lactic acid and trimethylene carbonate. These polymers have a long history of safe medical use and they degrade in vivo by hydrolysis into alpha-hydroxy acids that are metabolised by the body. Based on in-vitro testing, the **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** implants retain minimum of 70 % of their initial strength up to 9 weeks after implantation and gradually lose their strength thereafter. Bioresorption takes place within two to four years.

The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** implants are offered in different sizes and shapes and they are designed to be used with the **INION CPS™** instrumentation. The number of implanted plates and screws depends on the type of operation and the estimated maximum amount of implantable material to which patients can be exposed is 10 grams. The amount of implanted material constituents in this case is: L-lactic acid: 8,17 g; D-lactic acid: 1,07 g and Trimethylenecarbonate: 0,78g. This maximum exposure calculation includes 1 x **Inion CPS™** 14 x 14 hole mesh plate (4,7g) fixed with 30 x **Inion CPS™** 1,5 x 6 mm screws (á 0,02 g) and 1 x **Inion CPS™** 7 x7 hole mesh plate (1,2 g) + 10 x **Inion CPS™** plates (á 0,2-0,3 g) fixed with 74 x **Inion CPS™** screws (á 0,02-0,03 g), implanted in a single patient.

The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** implants are sterile, non-pyrogenic and non-collagenous.

INDICATIONS

Indications (CE mark):

The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** is intended for use in fractures, osteotomies and reconstructive procedures of the cranium, mid-face, orbital cavity, maxilla (e.g. LeFort I, II, III) and mandible in conjunction with appropriate maxillomandibular fixation.

Indications (USA):

- A. General indications: The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** is intended for use in trauma and reconstructive procedures in the craniofacial skeleton, mid-face, maxilla and mandible (in conjunction with appropriate maxillomandibular fixation). The **INION CPS™** Orbital Plates are intended for use in trauma and reconstructive procedures of the orbital cavity as part of the **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM**.
- B. Specific indications:
- Fractures of the cranium, mid-face, maxilla and mandible
 - Infant craniofacial surgery (i.e. craniosynostosis, congenital malformations)
 - LeFort (I, II, III) osteotomies
 - Pediatric reconstructive procedures
 - Orthognathic or reconstructive procedures of the cranium, mid-face, maxilla or mandible
 - Craniotomy flap fixation
- C. Additional indications: The **INION CPS™ 2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** meshes and screws (ref numbers PLT-1032, PLT-1033, PLT-1034, PLT-1035, SCR-1224, SCR-1225, SCR-1206, SCR-1207, SCR-1208, SCR-1290, SCR-1291, SCR-1292, SCR-1293, SCR-1294, SCR-1297, SCR-1298, SCR-1299, SCR-1300, SCR-1301, SCR-1209) are intended to maintain the relative position of bone grafts or bone graft substitutes in reconstructive procedures involving:
- Iliac crest harvest sites.

CONTRAINDICATIONS

- Mandibular tumor resection.
- Mandible without appropriate maxillomandibular fixation.
- Insufficient quality or quantity of bone.
- Active or potential infections.
- Patient conditions including limited blood supply, chronic disease which causes insufficient quality of bone, and where patient co-operation cannot be guaranteed (e.g. alcoholism, drug abuse).

TARGET POPULATION (CE mark):

Adult patients undergoing trauma or reconstruction procedures in the craniomaxillofacial area who have sufficient bone quality and quantity.

TARGET POPULATION (USA):

Patients undergoing trauma or reconstruction procedures in the craniomaxillofacial area who have sufficient bone quality and quantity.

WARNINGS

- Do not use the implant after expiration date printed on the label.
- Do not re-sterilize the implant.
- Do not re-use the implant.
- Do not use devices that are contained in opened or damaged packages.
- Use only if sterilization indicator shows exposed color red.
- Do not use if the temperature indicator dot is black.

- Remove the implant from the package using aseptic methods to retain implant sterility.
- Do not use for unintended applications! Proper function (i.e. effectivity and safety) of these implants cannot be guaranteed in case of off-label use.

PRECAUTIONS

- Only professional surgeons use the devices. The surgeon should be familiar with the devices, the method of application and the surgical procedure prior to performing the surgery.
- Instruments are available to aid accurate implantation of the **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM**. Surgical instruments are only to be used for their intended purpose. All instruments are to be regularly inspected for wear and damage. Use only the **INION CPS™** instruments for implant insertion.
- Use only **INION CPS™** screws to secure **INION CPS™** meshes and plates. Do not use **INION CPS™** screws to secure plates and meshes other than Inion's (e.g., metal plates).
- The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** devices provide fixation and are not intended to replace healthy bone or withstand the stress of full load bearing.
- Good alignment/reduction of the fracture/osteotomy is essential for optimal healing.
- Repeated bending can weaken the plate/mesh, especially once cooled. Warm bending can be carried out up to three times. Do not use plate(s)/mesh(es) that show signs of damage.
- Do not activate the screws by heating in the waterbath.
- Use only gentle pressure when placing screws. Do not overtighten.
- Using orbital plates without screws may lead to implant protrusion/gaze restriction/fibrous tissue formation.
- Meticulous hemostasis and complete primary soft tissue closure over the implant are essential.

RESIDUAL RISKS AND UNDESIRABLE SIDE-EFFECTS

Residual risks and undesirable side-effects with potential harm to the patient are similar to those with any method of internal fixation:

- Implantation of foreign materials can result in an inflammatory response or allergic reaction. Transient local fluid accumulation may occur in sterile circumstances, for example sterile abscess.
- Incorrect selection, placement, positioning or fixation of the implant can cause subsequent undesirable results, e.g. fixation failure, mechanical irritation.
- Fixation can break or loose as a result of early stress, activity or load bearing.
- Thin soft tissue coverage over the implant may increase the risk of complications.
- Implant may be covered with fibrous tissue capsule when implanted.
- Infections attributed to the surgical procedure.
- Neurovascular injuries that can occur due to surgical trauma.

Residual risks related to the insertion of the implants that may cause inconvenience for the user:

- Incorrect selection, placement, positioning, fixation or bending of the implant can cause breakage of the implant or instrument.
- Surgical instruments are subject to wear with normal usage and may break.

SPECIAL PATIENT POPULATIONS

- Patients with diabetes, bruxism, and temporomandibular joint malfunction or malformation may have an increased risk of complications.
- Patients with symptoms of sleep apnea have increased risk for complications following fixation of osteotomies in the maxilla and mandible.

SURGICAL INSTRUCTIONS

Implant selection

- The 1.5 System is recommended for use in cranial and mid-face fixation. Ø 2.0 mm screws are intended to be emergency screws for Ø 1.5 mm screws.
- The 2.0 System is recommended for use in mid-face and orthognathic maxillary fixation and mandibular trauma. The 2.0 System is not recommended for genioplasty. Ø 2.5 mm screws are intended to be emergency screws for Ø 2.0 mm screws.
- The 2.5 System is recommended for use in mandibular fixation and the 2.5 screws especially for genioplasty. Ø 2.8 mm screws are intended to be emergency screws for Ø 2.5 mm screws.
- The 2.8 BSSO screws are recommended for use in bicortical screw fixation of bilateral sagittal split osteotomies. Ø 3.1 mm screws are intended to be emergency screws for Ø 2.8 mm screws.

Preparation of the plates and meshes

The plates and meshes should be contoured to the bone according to these instructions (see also the instructions for use of the **INION THERMO™** waterbath):

- Place a sterile **INION THERMO™** drape to the **INION THERMO™** waterbath. DO NOT switch on the waterbath until filled with liquid!
- Pour sterile water or saline (min 0.4 litres max 0.7 litres) into the waterbath.
- Switch the waterbath on. Wait until the liquid heats up to the setting temperature 55°C / 131°F, this takes about 20 minutes.
- Choose appropriate **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** plate/mesh for the indication.
- Immerse the plate/mesh in the waterbath for at least a minute. After waterbath treatment: The plate/mesh is most malleable for 10-15 seconds, can be shaped as desired and may be cut to the desired size with standard surgical scissors.
- Maximum total immersion time of the plate/mesh in the waterbath is 30 minutes. Prolonged immersion may weaken the plate/mesh.
- Place the plate/mesh against bone at its intended location. Contour the plate/mesh to desired shape using fingers and/or the **INION CPS™** plate bending pliers. Take care to ensure that the plate/mesh is accurately adapted to the bone contours!
- If contouring is not satisfactory, the waterbath treatment can be repeated. If only part of the plate/mesh needs further adaptation dip only that portion of the plate/mesh.
- Repeated bending can weaken the plate, especially once cooled. Warm bending can be carried out up to three times. DO NOT use plate(s)/mesh(es) that show signs of damage!
- Take care not to perforate the **INION THERMO™** drape with surgical instruments.
- DO NOT activate the screws by heating in the waterbath!

Fastening of the plates and meshes manually with screws

- Choose appropriate **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** screws for fixation of the plate(s)/mesh(es). If needed, the screws can be cut to desired length with standard surgical scissors.
- DRILL & MANUAL TAP METHOD:

- Using the appropriate (corresponding to the screw diameter) **INION CPS™** bone drill bit attached to a slow speed drill (maximum speed of 2000 rpm), drill a suitable hole for the screw into the bone. Use irrigation during drilling. Always drill the holes for the screws through the holes of the plate/mesh to ensure correct location. Take care to hold the bone drill parallel to the long axis of the drill hole while drilling to prevent drill-wobble induced enlargement of the hole.
- Tap the screw holes through the holes of the plate/mesh manually using the appropriate (corresponding to the screw diameter) **INION CPS™** bone tap before screw insertion. Hold the bone tap parallel to the long axis of the drill hole both during insertion and reversing of the bone tap. Only use gentle pressure when tapping to preserve the prepared threads. Care must be taken to tap to a sufficient depth to allow complete seating of the screw. The prepared hole should be rinsed prior to screw insertion to flush out bone debris.
- **SELF-TAPPING DRILL METHOD:**
 - Using the appropriate (corresponding to the screw diameter) **INION CPS™** self-tapping drill bit attached to a slow speed drill, drill a suitable hole for the screw into the bone. Use irrigation during drilling. Always drill the holes for the screws through the holes of the plate/mesh to ensure correct location. Take care to hold the bone drill parallel to the long axis of the drill hole while drilling to prevent drill-wobble induced enlargement of the hole. The drilling should be stopped immediately before the drill stop contacts the plate or bone, and the drill bit should be reversed out carefully without stripping the prepared threads. Care must be taken to drill/tap to a sufficient depth to allow complete seating of the screw. The prepared hole should be rinsed prior to screw insertion to flush out bone debris.
- **SELF-TAPPING METHOD:**
 - If the bone is so thin or fragile, that tapping could cause breakage of the bone (e.g. maxilla), screws may be inserted in a self-tapping fashion. First use the appropriate (corresponding to the screw diameter) **INION CPS™** bone drill bit attached to a slow speed drill (maximum speed of 2000 rpm), drill a suitable hole for the screw into the bone. Use irrigation during drilling. Always drill the holes for the screws through the holes of the plate/mesh to ensure correct location. Take care to hold the bone drill parallel to the long axis of the drill hole while drilling to prevent drill-wobble induced enlargement of the hole. Take special care to only use gentle pressure when placing screws. Stop screwing once the screw head is seated in the plate hole. **DO NOT** overtighten!
- **SCREW LOADING & PLACEMENT:**
 - Align the appropriate **INION CPS™** screw driver with the screw head. Turn the screw driver until the screw rotates. Carefully but firmly push the screw driver into the screw head and lift the screw out rotating it slightly. Take care to ensure that the screw driver is aligned properly with the slot before pushing it into the screw, since there is a risk of damaging the plastic screw head.
 - Hold the driver and the screw parallel to the long axis of the drill hole and while keeping the plate/mesh in position, insert until the screw head is seated in the plated hole. **DO NOT** overtighten!
 - Screws can be used alone for the fixation of bone grafts and osteotomies. Hold the driver and the screw parallel to the long axis of the drill hole and while keeping the graft/bone fragments in position, insert the screw thread fully into the drill hole. **DO NOT** overtighten!
 - If the screw feels loose, or the screw head breaks during insertion: drill through the screw with the original size drill bit, tap the drill hole with one size larger bone tap and insert a corresponding larger diameter screw (i.e. emergency screw) instead.

POST-OPERATIVE REMINDERS

- As with any surgical procedure, careful postoperative management is important for optimal healing.
- Use appropriate maxillomandibular fixation especially with mandibular fractures and bimaxillary osteotomies (e.g. guiding elastics) during bone healing.
- Provide the patient with detailed instructions for postoperative care (e.g. hygiene maintenance, soft non chewing diet during bone healing, etc.).
- The patient should be warned that the fixation given with this device can break or loosen as a result of early stress, activity or load bearing.
- Antibiotic therapy at the discretion of the clinician.
- The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** is MR safe.
- Radiographs can be taken to check the alignment/reduction after fixation and to evaluate healing.
- Any serious incident that has occurred in relation to the **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** should be reported to the manufacturer and the competent authority of the EU Member State in which the user is established.

STERILITY

The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** implants have been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Do not use if the sterile package is opened or damaged. For single use only. Do not re-sterilize or re-use! Reprocessing may affect the properties of the implant (e.g. material degradation, implant deformation) and the safety and efficacy cannot be guaranteed. Do not use implant beyond the expiration date on the label.

STORAGE

Store at room temperature (15 to 30°C / 59 to 86°F). Protect from moisture. Product should not exceed a maximum temperature defined on the indicator. Do not use if the indicator dot is black.

DISPOSAL

Implants having been in contact with patient shall be disposed in accordance with hospital-specific instructions for biohazardous waste. Expired, unused implants shall be disposed in a controlled, secure manner to prevent misuse.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance (SSCP) of the **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** is available in the European database on medical devices (Eudamed). It is linked to the Basic UDI-DI of the product.

Eudamed public website: Not available yet

The **INION CPS™ 1.5/2.0/2.5 BIOABSORBABLE FIXATION SYSTEM** Basic UDI-DI: 06438408INIONCPSSYSTEMSD

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a licensed physician.

MANUFACTURER

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Finland
Tel. +358 10 830 6600
Fax +358 10 830 6601
info@inion.com
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

Medical device

Manufacturer

Date of manufacture

Keep dry

Do not re-sterilize

Use-by date

Double sterile barrier system

Unique device identifier

Sterilized using irradiation

Lot number

Catalogue number

Caution

Do not re-use

Storage temperature limits

Do not use if package is damaged and consult instructions for use

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İNION Oy
Inion® Inion Oy'un tescilli ticari markasıdır.