ENGLISH

**INION BIORESTORE™ (INION CPS BIORESTORE™ / INION GTR BIORESTORE™ / INION OTPS BIORESTORE™)
INION BIORESTORE™ PLUS**

INSTRUCTIONS FOR USE

DESCRIPTION

The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** is a bioactive, absorbable, synthetic bone graft substitute intended for filling, restoration, and/or augmentation of bony voids or gaps by new bone formation, which is achieved by the osteoconductive and osteostimulative properties of the porous bioactive material. The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** is used either separately or in conjunction with autogenous or allograft bone.

Bioactive materials are those materials that elicit a specific biological response at the interface of the material that results in the formation of a bond between the tissues and the material. Osteostimulation is a "...property of some bioactive materials to enhance, actively stimulate both the proliferation and differentiation of progenitor cells (e.g., mesenchymal stem cells)..."^{1,2}

¹ Society for Biomaterials presentation, 24th Annual Meeting, 1998:511-518

² J Appl Biomat 1992, 3:123-129

The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** implants are made of absorbable bioactive glass. Bioactive glasses have a long history of safe medical use and undergo a time-dependent, kinetic modification of the surface that occurs when implanted in living tissue. Specifically, the surface reaction results in the formation of a calcium phosphate layer that is substantially equivalent in composition and structure to the hydroxyapatite found in bone mineral. This apatite layer provides scaffolding onto which the patient's new bone will grow allowing complete repair of the defect. Based on preclinical testing and clinical data of **INION BIORESTORE™ / INION BIORESTORE™ PLUS**, most of the material remodels into bone and degrades *in vivo* in six months. The device lifetime is six months.

The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** consists of different size/shape porous block-shaped implants and morsels made of absorbable bioactive glass fibers. All products have the same effect on bone healing and bone formation as all the products are made of the same raw material composition and have an open porous structure. The block-shaped implants can be cut and shaped to desired size and shape according to the defect. The **INION BIORESTORE™ PLUS** block-shaped implants have slightly higher density for improved handling and shaping properties. Morsels are provided in nominal diameter of 3 mm. **INION CPS/OTPS/GTR BIORESTORE™** are additional trade names of **INION BIORESTORE™**.

Based on the maximum allowed volume of grafting material (30 cm³), considering the porosity of the implant, the maximum amount of implantable material to which a patient can be exposed is 39 grams. The amount of implanted material constituents in this case is: SiO₂ 19.7 g, K₂O 6.3 g, CaO 5.3 g, Na₂O 4.3 g, P₂O₅ 1.2 g, MgO 1.2 g, B₂O₃ 0.5 g, and TiO₂ 0.1 g.

The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** blocks may be cut or shaped to the desired three-dimensional shape with surgical instruments to fit the anatomy of surgical site. The **INION BIORESTORE™** morsels can be delivered into the defect either indirectly from the tube or directly from the syringe.

The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** products do not require specific instrumentation for their use.

The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** implants are sterile, non-pyrogenic and non-collagenous.

INDICATIONS

Indications (CE mark):

The **INION BIORESTORE™ (INION CPS BIORESTORE™, INION GTR BIORESTORE™, and INION OTPS BIORESTORE™)** morsels are indicated to be gently pressed into bony voids or gaps of the skeletal system, to fill, restore, and/or augment mandibular or maxillofacial bone defects and defects in the extremities, spine and pelvis. The osseous defects may be surgically created or resulting from traumatic injury to the bone.

The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** block-shaped implants are indicated to be gently pressed into bony voids or gaps of the skeletal system, to fill and/or restore defects in the extremities. The osseous defects may be surgically created or resulting from traumatic injury to the bone.

CONTRAINDICATIONS

- Defects that are intrinsic to the stability of the bony structure.
- Active or potential infection.
- Patient conditions including:
 1. Use of medication known to affect the skeleton (e.g., chronic glucocorticoid usage >10 mg/day for the previous 3 months). Estrogen replacement therapy is allowed.
 2. Need for chronic anticoagulant therapy (e.g., heparin). Prophylactic use of Coumadin or aspirin postoperatively is allowed.

3. A systemic metabolic disorder known to adversely affect bone healing and mineralization (e.g., poorly controlled insulin-dependent diabetes, renal osteodystrophy, Paget's disease), other than primary osteoporosis.
4. Any existing condition or disease that will interfere with good soft tissue and bone healing.
5. Defects requiring graft material to exceed the volume of 30 cm³.
6. Limited blood supply.
7. When patient cooperation cannot be guaranteed (e.g., alcoholism, drug abuse).
8. Patients with known allergy to the implant constituents or its degradation products.

TARGET GROUPS AND SPECIAL PATIENT POPULATIONS

The target patient population of the **INION BIORESTORE™ / INION BIORESTORE™ PLUS** consists of adult patients requiring bone grafting surgery and who are not contraindicated. There are no data from the use of the **INION BIORESTORE™** in pregnant women.

WARNINGS

- Do not use the implant after expiration date printed on the label.
- Do not re-sterilize the implant.
- Do not re-use the implant.
- Do not use devices that are contained in opened or damaged packages.
- Use only if sterilization indicator shows exposed color red.
- Remove the implant from the package using aseptic methods to retain implant sterility.
- Do not use for unintended applications! Proper function (i.e. effectivity and safety) of these implants cannot be guaranteed in case of off-label use.
- Do not use marrow in the syringe.
- Do not compress the material too strongly to prevent breakage of the material and its porous structure.
- Do not overfill the defect.
- Do not compromise blood supply to the defect area or blot away the blood/moisture in the placed graft material.
- Do not leave defect open.

PRECAUTIONS

- Only professional surgeons use the devices. The surgeon should be familiar with bone grafting and internal/external fixation techniques, the devices, the method of application and the surgical procedure prior to performing the surgery.
- The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** devices are not intended to withstand the stress of immediate load bearing prior to soft and hard tissue ingrowth without appropriate additional fixation. In cases of fracture fixation or load bearing applications, standard internal or external stabilization techniques must be followed to obtain rigid stabilization in all planes. Instrumentation used in conjunction with the **INION BIORESTORE™ / INION BIORESTORE™ PLUS** must gain purchase in the host bone. The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** may not be used to gain immediate screw purchase or to stabilize screw placement. Mechanical loading should not take place for at least 5-6 months after grafting.
- Block preparation: Take care to ensure that any shaped device surfaces are smooth and free from excessive loose particles.
- Thin soft tissue coverage over the implant, or use in areas where the graft cannot be adequately contained (to prevent motion and migration of the material) may increase the risk of complications.
- Syringe: Always protect the mixture from contamination by keeping the cap of the syringe on and keeping the syringe inside the sterile field until the mixture has been delivered into the defect.

RESIDUAL RISKS AND UNDESIRABLE SIDE-EFFECTS

Residual risks and undesirable side-effects with potential harm to the patient are similar as to be expected of autogenous bone grafting procedures:

- Implantation of foreign materials can result in an inflammatory response or allergic reaction.
- Incorrect selection, placement or positioning of the implant can cause subsequent undesirable results (e.g. delayed bone formation, mechanical irritation), or breakage of implants.
- Implant can break or loosen as a result of early stress, activity or load bearing.
- Infections attributed to the surgical procedure, such as superficial wound infection, deep wound infection, and deep wound infection with osteomyelitis.
- Neurovascular injuries that can occur due to surgical trauma.
- Thin soft tissue coverage over the implant, or use in areas where the graft cannot be adequately contained (to prevent motion and migration of the material) may increase the risk of complications.
- Other complications that may arise as a result of surgery may include: delayed bone formation, failure of restoration, loss of bone graft, graft protrusion and/or dislodgement, pain, swelling, and general complications that may arise from anesthesia and/or surgery.

SURGICAL INSTRUCTIONS

Implant selection

All **INION BIORESTORE™ / INION BIORESTORE™ PLUS** product variants can be used according to the indications for use based on surgeon's selection as suitable to the particular site. Different packaging volumes of morsels are available, packaged in a tube or in a syringe. Different sizes and shapes of block-shaped implants are available and they are intended to be shaped to fit the particular defect by the surgeon.

1. Surgical considerations and reminders

- Prophylactic preoperative antibiotic treatment is recommended.
- Use proper local, regional or general anaesthesia.

- Maintain sterile field throughout the procedure.
- Proper exposure using standard surgical procedure.
- Thoroughly prepare the surgical site preserving the neurovascular structures by careful dissection.
- Make every attempt to obtain complete tension-free primary soft tissue closure over the material.

2. Preparation

- Choose appropriate **INION BIORESTORE™ / INION BIORESTORE™ PLUS** block(s) or morsels for the indication. Preoperative radiographic evaluation of the defect site is essential to accurately assess the extent of the defect and to aid in the selection and placement of the graft material.
- SHAPING OF BLOCKS:
 - The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** blocks may be cut or shaped to the desired three dimensional shape with surgical instruments (e.g., surgical knives, scalpels or scissors). However, before shaping the **INION BIORESTORE™ / INION BIORESTORE™ PLUS**, the implant should be moistened with either sterile saline or water. Furthermore, the shaped implant should always be rinsed with sterile saline or water after shaping and before any contact with the patient's tissues to remove excessive loose particles from the surfaces of the device. Take care to ensure that any shaped device surfaces are smooth and free from excessive loose particles.
- WETTING:
 - The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** should be wet with patient's blood or marrow, or sterile saline or water before application. The osseous regeneration may occur more rapidly by using the patient's osteogenic blood or marrow obtained by scoring the surface or drilling into the cortical bone at the surgical site. The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** forms a cohesive mass when mixed with patient's blood. Place the desired fluid in a sterile basin or other suitable container, or inject blood directly into the syringe containing **INION BIORESTORE™** morsels. Place the material in contact with the fluid to soak up the fluid. The porous structure of the material will allow the fluid to permeate into the implant by capillary action.
- SYRINGE PREPARATION:
 - Pull the piston of the syringe back to be fully extended. Turn the syringe upright (cap upwards), gently tap the container to release the morsels from the cap, and open the cap.
 - Add the patient's blood carefully into the syringe, and close the cap. Do not use marrow in the syringe.
 - Turn the syringe slowly around to wet the morsels completely. If needed, move the piston slowly back and forth to enhance moistening, but be cautious not to pressurize the container excessively. Allow to set until the morsels have soaked up the blood and the mixture has thickened as desired. The mixture can be delivered immediately after the morsels have soaked up the blood, but if the mixture is allowed to set the mixture will become thickened which may be desirable depending on the size and shape of the defect. The maximum recommended immersion time is 60 minutes. Always protect the mixture from contamination by keeping the cap of the syringe on and keeping the syringe inside the sterile field until the mixture has been delivered into the defect.
 - Open the cap before implantation.

3. Implantation

- SITE PREPARATION:
 - Once the defect has been exposed, thorough debridement of the osseous defect of all necrotic and soft tissue should be carried out. Eliminate all granulation tissue or cysts at the defect site. The walls of the defect should be decorticated.
 - Irrigate the defect with sterile saline or water and evacuate the excess.
 - Some bleeding should be observed originating from the host bone to indicate viability. Intra-marrow penetration is useful to ensure some bleeding from the host bone, which provides a supply of osteogenic material to aid in bone regeneration.
- PROPER PLACEMENT AND CONTAINMENT OF THE DEVICE:
 - Place the prepared **INION BIORESTORE™ / INION BIORESTORE™ PLUS** bone graft substitute into the intended graft site. The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** should fill the defect and contact viable bone as much and as tightly as possible. Regeneration will occur best when the material is in tight contact with the bony walls of the defect, and blood and blood vessels can infiltrate the graft material. However, when placing the bone graft substitute into the defect, **DO NOT** compress the material too strongly to prevent breakage of the material and its porous structure!
 - The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** should be used in areas where the graft can be adequately contained to prevent motion and migration of the material.
 - The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** should fill the bony defect to the highest level of the defect. However, **DO NOT** overfill the defect! **If the selected implant is too large for the defect, the implant should be shaped to the desired shape according to the above described instructions.**
- SITE CLOSURE:
 - After placement of the **INION BIORESTORE™ / INION BIORESTORE™ PLUS**, complete tension-free primary soft tissue closure over the grafted area.
 - **DO NOT** compromise blood supply to the defect area or blot away the blood/moisture in the placed graft material!
 - **DO NOT** leave defect open!
 - Degradable or non-degradable membranes or graft containment devices may also be used to achieve closure. The use of a membrane or graft containment device is recommended in cases where the defect is large or limited bony retention is present.

These instructions are intended as guideline for the use of the **INION BIORESTORE™ / INION BIORESTORE™ PLUS** implants as a part of established techniques. They are not intended to replace or change standard grafting techniques.

POST-OPERATIVE REMINDERS AND INFORMATION OF CLINICAL BENEFITS

- As with any surgical procedure, careful postoperative management is important for optimal healing. Postoperative patient management should follow the same regimen as similar cases utilizing autogenous bone grafting. Standard postoperative practice should be followed.
- Use appropriate additional immobilization.
- Provide the patient with detailed instructions for postoperative care (e.g. regarding immobilization and hygiene maintenance).

- The patient should be warned that the implants can break or loosen as a result of early stress, activity or load bearing to ensure reduced loading to prevent collapse and deformity.
- After treatment of mandibular or maxillary bone defects, the patient should be cautioned against chewing on or brushing the area for at least 1 to 2 weeks or as long as the clinician deems necessary. Antibacterial rinses should be used during this period.
- Antibiotic and/or analgesic therapy at the discretion of the clinician.
- Any serious incident that has occurred in relation to the **INION BIORESTORE™ / INION BIORESTORE™ PLUS** should be reported to the manufacturer and the competent authority of the EU Member State in which the user is established.
- The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** provides several benefits for bone grafting surgery:
 - Bone healing (formation/fusion/restoration/augmentation) enhanced with osteostimulative, bioactive material:
The bioactive, osteostimulative properties of **INION BIORESTORE™ / INION BIORESTORE™ PLUS** accelerate the regeneration rate of new tissue by actively stimulating osteoblast proliferation and differentiation, while providing a scaffold to support new bone formation. The formation of fibrous tissue is avoided with the use of bioactive glass materials. These materials allow the bone to bond directly on the implant surface. The bioactive, osteostimulative properties have been demonstrated in *in vitro* degradation and cell culture tests.
 - Controlled absorption and remodelling to bone:
The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** gradually absorbs and remodels into new bone in an efficient manner due to its porous fiber structure. This open structure enables flow of fluid into the implants, and allows bone formation and vascularization throughout the implant.
 - Avoidance of donor site surgery/morbidity/complications; reduced total operative time, hospitalization, and anesthesia:
INION BIORESTORE™ / INION BIORESTORE™ PLUS enables a single procedure in the recipient site only, in reduced operative time.
 - Antimicrobial, synthetic, sterile:
Based on laboratory testing, the **INION BIORESTORE™ / INION BIORESTORE™ PLUS** has antimicrobial properties, which may inhibit the occurrence of bacterial infection at the surgical site. The synthetic origin of the **INION BIORESTORE™ / INION BIORESTORE™ PLUS** excludes risks of disease transmission. The overall adverse event risk is reduced.
 - No imaging interference:
No artefact is created on computer tomography (CT) and magnetic resonance imaging (MRI) by the material.
 - The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** is MR safe. The **INION BIORESTORE™ / INION BIORESTORE™ PLUS** bone graft substitute is entirely made of non-conductive and non-ferrous material components. Measured levels of metallic impurities are negligible.

STERILITY

The **INION BIORESTORE™** has been sterilized with ionizing irradiation. Use immediately after opening the sterile seal. Do not use if the sterile package is opened or damaged. For single use only. Do not re-sterilize or re-use! Reprocessing may affect the properties of the implant (e.g. material degradation, implant deformation) and the safety and efficacy cannot be guaranteed. Do not use implant beyond the expiration date on the label.

STORAGE

Protect from moisture.

DISPOSAL

Implants having been in contact with patient shall be disposed in accordance with hospital-specific instructions for biohazardous waste. Expired, unused implants and remainders shall be disposed in a controlled, secure manner to prevent misuse.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The summary of safety and clinical performance (SSCP) of the **INION BIORESTORE™ / INION BIORESTORE™ PLUS** is available in the European database on medical devices (Eudamed). It is linked to the Basic UDI-DI of the product.

Eudamed public website: Not available yet.

The **INION BIORESTORE™** Basic UDI-DI: 6438408INIONBIORESTORERQ.

MANUFACTURER

INION Oy
Lääkärintäti 2
33520 Tampere
Finland
Tel. +358 10 830 6600
info@inion.com
www.inion.com

CE 2797

Inion® is a registered trademark of Inion Oy.

SYMBOLS USED IN LABELING

Medical device

Manufacturer

Date of manufacture

Keep dry

Do not re-sterilize

Use-by date

Double entry package with two sterile barrier systems

Single sterile barrier system (syringe)

Unique Device Identifier

Sterilized using irradiation

Lot number

Catalogue number

Caution, consult accompanying documents

Do not re-use

Do not use if package is damaged and consult instructions for use

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¹ Präsentation der Society for Biomaterials, 24th Annual Meeting, 1998:511-518

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Inion® on Inion Oy:n rekisteröity tavaramerkki.