ENGLISH

## **INSTRUCTIONS FOR CLEANING, DISINFECTION AND STERILIZATION OF REUSABLE INION® INSTRUMENTATION**

Cleaning, disinfection and sterilization are required before the first use of non-sterile **INION® INSTRUMENTATION** and before each re-use of all reusable **INION® INSTRUMENTATION**. Effective cleaning and disinfection are a vital requirement for an effective sterilization of the **INION® INSTRUMENTATION**. Cleaning, disinfection and sterilization of reusable **INION® INSTRUMENTATION** must be conducted according to the procedures outlined in this document.

### **HOW SUPPLIED**

**INION® INSTRUMENTATION** is supplied as nonsterile, and must be cleaned and sterilized prior to use according to the procedures outlined in this document.

### **WARNINGS**

- Do not expose the **INION® INSTRUMENTATION** to temperatures higher than 140°C (284°F).
- Do not use high **acidic or high alkaline solutions**.<sup>1)</sup>
- Do not use saline solutions.
- Do not use any instrument as a lever arm.
- Use the instrument only for its intended use.
- Do not change limb position while instrument is in the joint.
- Excess load directed to the instruments may cause breakage of the instrument tip.
- Careless use of instruments with sharp edges may result in soft tissue damage.

### **INSTRUCTIONS: POINT OF USE**

- Do not allow blood or soil to dry on the instruments.
- Remove excess soil with disposable, non-shedding wipe.
- Flush cannulated devices with cool water.
- It is recommended that instruments are reprocessed as soon as is reasonably practical following use.

### **MANUAL CLEANING/DISINFECTION**

#### **PRE-CLEANING:**

Equipment: Soaking bath, detergent, brush, running water/purified water, ultrasonic washer unit, DI water, syringes min. volume 20 ml.

Ensure the compatibility of the cleaning agent with the **INION® INSTRUMENTATION**.

Follow the cleaner or detergent manufacturer's instructions for use for correct exposure time, temperature, water quality and concentration.

1. Instruments with removable parts must be disassembled prior to cleaning.
2. Immerse the instrument in a non-foaming, neutral enzymatic cleaner or neutral detergent (7.0-9.0 pH) solution and allow to soak for a minimum of 10 minutes. Alkaline agents with a pH of 12 or less may be used to clean stainless steel and polymer instruments in countries where required by law or local ordinance; or where prion diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeld-Jakob Disease (CJD) are a concern.<sup>1)</sup> This treatment does not have adverse effect to stainless steel instruments. Aluminium and plastic parts may suffer visual deterioration.
3. Clean the instrument with a soft-bristled brush and common hospital cleaning tools for a minimum of five minutes. Pay extra attention to challenging design features and areas, where tissue residue accumulates the most easily, and for cannulated instruments, to the cavity/area, which overlaps with a potential rinsing adapter. Move and/or retract all moveable device features and use brush to clean hinges and mating surfaces. Scrub lumens and holes with a brush. Clean long, narrow cannulations meticulously with a brush or by injecting water through the cannula.
4. Rinse under clean running water for min. 2 minutes. Ensure that running water passes through cannulations and that blind holes are repeatedly filled and emptied. Movable features should be actuated during rinsing. It is critical that alkaline cleaning agents are completely and thoroughly neutralized and rinsed from the devices or degradation may occur that limits the device life.
5. Sonicate the instrument in warm enzymatic detergent in its fully opened position for a minimum of 10 minutes.
6. Rinse the instrument thoroughly with warm deionized/purified water for a minimum of 3 minutes. Ensure that running water passes through cannulations and that blind holes are repeatedly filled and emptied. Rinse all lumens of the instruments once by application of a single-use syringe. For cannulated instruments with handle, connect syringe directly to the Luer cavity at the handle for enhanced water flow. For other cannulated instruments, use suitable adapter or aim for concentrated water flow to the cannula. Movable features should be actuated during rinsing.

<sup>1)</sup>The capacity of a detergent to inactivate prions is most likely with pH-values > 10 and a soaking time of more than 10 minutes under elevated, but not protein-fixing temperatures (for example 55°C). However, claims that a product can inactivate prions must be evidence-based. Hygiene Requirements for the Reprocessing of Medical Devices, Robert Koch Institute and BfArM. Bundesgesundheitsbl 2012 • 55:1244–1310.

#### **DISINFECTION:**

Equipment: Disinfection bath, disinfectant, syringes min. volume 20 ml, deionized water.

Ensure the compatibility of the disinfectant with the **INION® INSTRUMENTATION**.

Follow the disinfectant manufacturer's instructions for use for correct exposure time, temperature, water quality and concentration.

1. Soak the disassembled products for a minimum of 10 minutes in the disinfectant solution so that the products are sufficiently covered. Pay attention that there is no contact between the products. Sway movable parts several times during disinfection.
2. Rinse all lumens of the products at least five times at the beginning of the soaking time by application of a single-use syringe. For cannulated instruments with handle, connect syringe directly to the Luer cavity at the handle for enhanced water flow. For other cannulated instruments, use suitable adapter or aim for concentrated water flow to the cannula.
3. Remove the products of the disinfectant solution and post-rinse them at least five times by soaking (at least 1 min) with deionized water.
4. Rinse all lumens of the products at least five times by application of a single-use syringe. For cannulated instruments, see step 2 above.
5. Dry the instruments thoroughly by clean, soft, lint-free cloth. Move and/or retract all movable parts during drying. Clean compressed air may be used for drying of inaccessible features. If additional drying is required, place the instruments in a clean area or heat in an oven below 90°C (194°F).

The fundamental suitability of the **INION® INSTRUMENTATION** for an effective manual cleaning and disinfection was demonstrated by an independent, governmentally accredited (DAkkS, with ILAC Mutual Recognition Agreement) in accordance with DIN EN ISO/IEC 17025 and GLP certified test laboratory by application of an ultrasonic bath of the SONOREX series with 35 kHz (BANDELIN electronic, D-Berlin) for manual cleaning as well as the cleaning detergent Cidezime/Enzol and the disinfectant Cidex OPA (ASP, distributed by Johnson & Johnson, D-Norderstedt) considering the specified procedure (but application of worst case settings with respect to the instructions of the detergent manufacturer).

## AUTOMATED CLEANING/DISINFECTION

Automatic washer/disinfector using thermal disinfection is recommended for **INION® INSTRUMENTATION**. Follow the internal hospital instructions and recommendations of manufacturers of automatic washer/disinfector, especially in regard to operating time and temperature. Equipment: Automatic washer/disinfector, cannula brush, running water, cleaner, detergent, deionized or purified water.

Ensure the compatibility of the cleaning agent with the **INION® INSTRUMENTATION**.

- Do not expose to temperatures higher than 140°C (284°F).
- **INION® INSTRUMENTATION** is classified as critical according to the Ao classification as defined in ISO 15883 Washer-disinfectors.
- Use automatic washer/disinfector according to EN ISO/ANSI AAMI ST15883 and with fundamentally approved efficiency (for example CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration).
- Sterile / deionized / purified water is recommended for final rinsing.
  1. Conduct steps 1-6 above from section Manual Pre-cleaning.
  2. Active rinsing of instruments with cannula diameter < 4 mm: Instruments with a fixed handle – connect with Luer adapter (e.g. Miele E791) and potential extension (e.g. Miele E792). Instruments with AO or other coupling – connect with AO/Luer adapter and potential extension (e.g. Miele E792), or place instrument upright to a rinsing tube (e.g. Miele E442). Instruments with cannula diameter ≥ 4 mm: place upright directly to the rinsing shaft of washer/disinfector.
  3. Place the instrument in the automatic washer/disinfector so that the instrument cannulations, hinges, lumens and holes can drain and use a washing/disinfection cycle with the following minimum properties:

Run the instruments through a standard instrument cycle of the washer/ disinfector by selecting an approved program for thermal disinfection (**A<sub>0</sub> value ≥ 3000 or – in case of older devices - at least 5 min at 90°C / 194°F**), with a sufficient number of rinsing steps (**at least three** rinsing steps after cleaning / respectively neutralization, if applied) or conductance based rinsing control recommended in order to effectively prevent remnants of the detergents.

Post-rinsing only with sterile or low contaminated water (max. 10 germs/ml, max. 0.25 endotoxin units/ml), for example purified/highly purified water.

Run the automatic drying process of the washer/disinfector min. 30 minutes at 60°C / 140°F. If necessary, additional manual drying with lint-free cloth and lumen insufflation with filtered compressed air (oil-free, low contamination with microorganisms and particles).

**ALTERNATIVELY**, use a washing/disinfection cycle with the following minimum parameters:

Pre-wash	5 minutes	Cold drinking water, < 43°C (109°F)
Wash I	5 minutes	Neutral pH (7.0-9.0) enzymatic cleaner, warm drinking water (> 45°C / 113°F)
Wash II	5 minutes	Neutral pH (7.0-9.0) detergent wash, warm drinking water (> 45°C / 113°F)
Rinse	5 minutes	Warm DI or PURW (> 40°C / 104°F)
Thermal disinfection	5 minutes	90°C (194°F) <sup>2)</sup>
Drying	40 minutes	

The fundamental suitability of the **INION® INSTRUMENTATION** for an effective automated cleaning and disinfection was demonstrated by an independent, governmentally accredited (DAkkS, with ILAC Mutual Recognition Agreement) in accordance with DIN EN ISO/IEC 17025 and GLP certified test laboratory by application of the detergent Cidezime (ASP, distributed by Johnson & Johnson, D-Norderstedt) and an ultrasonic bath of the SONOREX series with 35 kHz (BANDELIN electronic, D-Berlin) for manual pre-cleaning as well as the WD G 7836 CD, Miele & Cie. GmbH & Co., D-Gütersloh (thermal disinfection), the program DES-VAR-TD, and the cleaning detergent Neodisher medizym (Dr. Weigert GmbH & Co. KG, D-Hamburg) considering to the specified procedure (but application of worst case settings with respect to the program DES-VAR-TD and the instructions of the detergent manufacturer).

## INSPECTION

- Instruments should be carefully inspected prior to sterilization in order to ensure that all visible contamination has been removed.

## REASSEMBLY

- If appropriate, reassemble instruments before sterilization.

## FUNCTIONAL TESTING

- Actuate the moving parts to ensure smooth operation throughout the intended range of motion.
- If the instrument is part of a larger assembly, check that devices can be assembled with mating components.
- Examine the instruments for signs of wear, corrosion, bluntness and non-functioning parts to ensure proper performance of the device.

## STERILE PACKAGING

- Packaging should ensure sterility of the instruments until opened for use at the sterile field.
- The **INION® INSTRUMENTATION** must be packaged separately into a single-use sterilization packaging suitable for steam sterilization or in the **INION®** sterilization tray designed for the purpose.

## STERILIZATION

- Steam sterilization with pre-vacuum procedure is recommended for the **INION® INSTRUMENTATION**.
- The **INION® INSTRUMENTATION** has to be cleaned and disinfected prior to sterilization.

The recommended steam sterilization parameters are as follows<sup>2)</sup>:

Sterilization method	Temperature	Exposure time	Configuration	Drying time
Pre-vacuum	121°C / 250°F	20 minutes	In sterilization pouch	20 min in pouch
Pre-vacuum	132°C / 270°F	3 minutes	In full tray with sterilization wrap	60 min in wrap tray cover off
Pre-vacuum	134°C / 273°F	3 minutes	In full tray with sterilization wrap	60 min in wrap tray cover off
US FDA: Pre-vacuum	132°C / 270°F	4 minutes	In full tray with sterilization wrap	60 min in wrap tray cover off
Germany: Pre-vacuum	134°C / 273°F	5 minutes	In full tray with sterilization wrap	60 min in wrap tray cover off

<sup>2)</sup> Providing the medical instrument is heat resistant, an autoclave cycle at 134°C for 18 minutes or more is generally effective for prion inactivation (Fichet et al. 2004: Novel methods for disinfection of prion-contaminated medical devices. Lancet 364, 521–6.)

## STORAGE

- Store sterilized instruments in a cool, dry place protected from dust and direct sunlight.
- Examine carefully prior to opening the package that the package integrity has not been compromised.

## DEVICE LIFETIME AND LIMITATIONS ON REPROCESSING

Repeated processing according to the instructions has a minimal effect on the **INION® INSTRUMENTATION**. Life cycle is normally determined by wear and damage due to use. Damaged **INION® INSTRUMENTATION** must not be further used and should be disposed of. These instructions have been validated by Inion Oy as being capable of preparing the **INION® INSTRUMENTATION** for re-use for 100 processing cycles. For additional processing cycles, the end user is responsible for ensuring proper functionality, durability and biocompatibility of the device by validation and/or periodically performed checks.

The end user is responsible for the sterility of the **INION® INSTRUMENTATION**. Any deviation from these instructions should be properly evaluated for effectiveness and potential adverse consequences. Follow the internal hospital instructions and recommendations by manufacturers of detergents and of automatic washers and sterilizers.

## MANUFACTURER

INION Oy  
Lääkärinkatu 2  
33520 Tampere  
Finland  
Tel. +358 10 830 6600  
info@inion.com  
www.inion.com

CE / CE 2797

Inion® is a registered trademark of Inion Oy.

## SYMBOLS USED IN LABELING

Medical device

Unique Device Identifier

Manufacturer

Non-sterile product

Consult instructions for use

Lot number

Catalogue number

Date of manufacture

Keep dry

Caution, consult accompanying documents

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Inion® Inion Oy'un tescilli markasıdır.