



You noticed me
because of my color



Visibly Different
Inion Trinion™ The first
colored biodegradable
meniscus screw

Other visible advantages of the Inion Trinion™ meniscus screw: It is flexible and durable, dual threaded for controlled tear closure, unique self-tapping tip for smooth and easy insertion, headless screw design to avoid adjacent tissue damage. Trinion™ gives optimal insertion control, remaining visible until it is fully inserted within the meniscus.

Inion Ltd: tel +358-3-2306 600 • Inion Inc: Toll-free telephone (866-Inion-US) • info@inion.com | www.inion.com

