


Inion FixOn™

Titanium Suture Anchor



Inion FixOn™ - Titanium Suture Anchor

The safe and reliable Inion FixOn™ Ti Suture Anchors consist of a titanium suture anchor attached with strong, non-absorbable braided sutures preloaded on a practical driver.

Anchors are available from Ø2.0 mm to Ø5.5 mm and thus provide a solution for a wide range of surgical operations. Inion FixOn™ Ti Suture Anchors are sterile and intended for single use only.



Benefits

- **Excellent fixation strength**
The Inion FixOn™ Ti Suture Anchors provide safe and reliable fixation for orthopedic surgical operations in knee, shoulder, elbow, ankle, foot, hand and wrist.
- **High-quality anchor material**
The suture anchors are made of a high-purity titanium alloy (Ti6AlV ELI).
- **Strong, clearly visible sutures (white/striped/blue)**
 - Suprabraid™ sutures are prepared from strong, ultrahigh molecular weight polyethylene (UHMWPE) that is easy to knot and handle.
 - Anchors are provided with #2 (5 metric) sutures, except for the Ø2.0 mm anchor which is provided with #2-0 (3 metric) suture.
- **Fast and easy insertion**
The Ø3.5, Ø4.5, and Ø5.5 mm suture anchors can be simply inserted with a “punch and drive” technique.



Indications (CE mark)

Inion FixOn™ Ti Suture Anchors are intended to reattach soft tissue to bone in orthopedic surgical operations (i.e. arthroscopic or open surgical procedures) in knee, shoulder, elbow, ankle, foot, hand and wrist where the anchor size is appropriate in the presence of appropriate postoperative immobilization throughout the healing period.


Ordering Information

Implants:

FOT-2001	Ø2.0 mm Suture Anchor with one suture
FOT-2801	Ø2.8 mm Suture Anchor with one suture, long shaft
FOT-3502	Ø3.5 mm Suture Anchor with two sutures
FOT-4502	Ø4.5 mm Suture Anchor with two sutures
FOT-4503	Ø4.5 mm Suture Anchor with three sutures
FOT-5502	Ø5.5 mm Suture Anchor with two sutures
FOT-5503	Ø5.5 mm Suture Anchor with three sutures

Instruments:

INS-9637	Drill bit Ø1.2 mm1
INS-9638	Awl
INS-9639	Guide
INS-9640	2nd level tray (INS-9640 can be accommodated in INS-9619 Inion® Shoulder instrument tray)



INION OY

Lääkärintäti 2
FI-33520 Tampere, FINLAND
tel: +358-10-830 6600
fax: +358-10-830 6601
email: info@inion.com
internet: www.inion.com

INION INC

2800 Glades Circle
Suite 138, Weston
FL 33327
USA
Toll-free tel: 866-INION-US
tel: 954-659-9224-fax: 954-659-7997