


Instrument Tray

Inion CompressOn™

INION



Inion CompressOn™ instrument tray and instrumentation are designed to be used with the Inion CompressOn™ bioabsorbable headless screw. The instruments are CE-marked medical devices, manufactured by Inion Oy. The conformity assessment was conducted by BSI The Netherlands (CE 2797).

Inion CompressOn™ Instrument Tray

The Inion CompressOn™ Instrument Tray is a compact and easy-to-use tray for various sizes of instrumentation. Inion CompressOn™ instrumentation is designed to be used with bioabsorbable Inion CompressOn™ screws. The instrumentation includes bone drill bits, drill taps, screwdriver shafts, screwdriver handles, drill sleeves, K-wires, and depth gauges.



All parts of the Inion CompressOn™ Instrument Set that are in direct patient contact are manufactured from stainless steel materials suitable for surgical instruments. Sterilization tray instrument holders are manufactured from medical-grade silicone or stainless steel.

The Inion CompressOn™ tray comes with plastic identification plates that can be easily changed depending on the size of the instrumentation needed.



Product Overview

The following items are available for the Inion CompressOn™ Instrument Tray IFS-9802:

Art. No	Screw size	Description	Art. No	Screw size	Description
IFS-2703	2.7	Drill bit Ø 2.0 mm, cannulated	IFS-9041	2.7 / 3.5 / 4.5 / 5.0	Depth gauge Ø 2.7/3.5 mm
IFS-2708	2.7	Drill tap Ø 2.7 mm, cannulated	IFS-9044	6.0	Depth gauge for 230 mm K-wire
IFS-3503	3.5	Drill bit Ø 2.7 mm, cannulated	IFS-9045	2.7 / 3.5 / 4.5 / 5.0	Depth gauge for 150 mm K-wire
IFS-3508	3.5	Drill tap Ø 3.5 mm, cannulated	IFS-9063	2.7	K-wire Ø 0.8 x 150 mm (10 pieces)
IFS-4503	4.5	Drill bit Ø 3.2 mm, cannulated, long	IFS-9064	3.5	K-wire Ø 1.25 x 150 mm (10 pieces)
IFS-4508	4.5	Drill tap Ø 4.5 mm, cannulated, long	IFS-9065	4.5 / 5.0	K-wire Ø 1.6 x 150 mm (10 pieces)
IFS-4513	4.5	Drill bit Ø 3.2 mm, cannulated, short	IFS-9066	6.0	K-wire Ø 2.0 x 230 mm (10 pieces)
IFS-4518	4.5	Drill tap Ø 4.5 mm, cannulated, short	INS-9093	2.7 / 3.5 / 4.5 / 5.0 / 6.0	Small handle, cannulated
IFS-5003	5.0	Drill bit Ø 3.7 mm, cannulated	INS-9120	4.5 / 5.0 / 6.0	Large handle, cannulated
IFS-5008	5.0	Drill tap Ø 5.0 mm, cannulated	INS-9121	4.5 / 5.0 / 6.0	Large handle, ratchet, cannulated
IFS-6003	6.0	Drill bit Ø 4.2 mm, cannulated	IFS-9061	3.5	K-wire Ø 1.25 x 230 mm (10 pieces)
IFS-6008	6.0	Drill tap Ø 6.0 mm, cannulated	IFS-9062	4.5 / 5.0	K-wire Ø 1.6 x 230 mm (10 pieces)
IFS-9003	2.7	Screwdriver shaft, T7, cannulated	IFS-9043	4.5 / 5.0	K-wire Ø 1.6 x 300 mm with depth gauge
IFS-9004	3.5	Screwdriver shaft, T8, cannulated	IFS-9067	3.5	K-wire Ø 1.25x150 mm, threaded
IFS-9005	4.5 / 5.0	Screwdriver shaft, T15, cannulated	IFS-9068	4.5 / 5.0	K-wire Ø 1.6x150 mm, threaded
IFS-9006	6.0	Screwdriver shaft, T25, cannulated	IFS-9069	6.0	K-wire Ø 2.0x230 mm, threaded
IFS-9022	2.7	Drill sleeve Ø 2.0-2.7 mm			
IFS-9023	3.5	Drill sleeve Ø 2.5-3.5 mm			
IFS-9024	4.5	Drill sleeve Ø 3.2-4.5 mm			
IFS-9025	5.0	Drill sleeve Ø 3.7/5.0 mm			
IFS-9026	6.0	Drill sleeve Ø 4.2/6.0 mm			

Please confirm up-to-date product availability with Inion Sales.



INION OY
Lääkärintäti 2
FI-33520 Tampere
FINLAND
tel: +358-10-830 6600
email: info@inion.com
internet: www.inion.com

INION INC
37 North Orange Avenue, Suite 900C
Orlando, FL 32801
USA
Toll-free tel: 866-INION-US
tel: +1 954-659-9224
email: inion.us@inion.com

Version 2: 9/2026