


SECURE THE K-WIRE IN PLACE

INION

Developed to address unintended backing out, our threaded K-wire incorporates a retention mechanism, helping to maintain **secure engagement** and consistent control throughout the surgical procedure. Compatible with our Inion CompressOn™ screws, it supports procedural efficiency, surgical confidence, and reliable fixation.



Threaded K-Wire

as part of the Inion CompressOn™ instruments

Prevents the K-wire from inadvertently detaching from the bone during surgery.

Threaded K-Wires product information

IFS-9067	K-wire Ø1.25x150 mm, threaded	Ø3.5 mm CompressOn™ screws
IFS-9068	K-wire Ø1.6x150 mm, threaded	Ø4.5 mm CompressOn™ screws Ø5.0 mm CompressOn™ screws
IFS-9069	K-wire Ø2.0x230 mm, threaded	Ø6.0 mm CompressOn™ screws

INION OY
Lääkärintäti 2
33520 Tampere

✉ info@inion.com

www.inion.com
v1 7/2026